Tokai Microcatheter

Camelian

* B M P 0 5 3 E *

Camelian
Support

CE 1639

INSTRUCTIONS FOR USE

English.....	2	Lietuviškai.....	44
български.....	5	Norsk.....	47
Hrvatski.....	8	Polski.....	50
Česky.....	11	Português.....	53
Dansk.....	14	Română.....	56
Nederlands.....	17	Русский.....	59
Eesti.....	20	Srpski.....	62
Suomi.....	23	Slovensky.....	65
Français.....	26	Slovenščina.....	68
Deutsch.....	29	Español.....	71
Ελληνικά.....	32	Svenska.....	74
Magyar.....	35	Türkçe.....	77
Italiano.....	38	Symbols.....	80
Latviski.....	41		

TOKAI MEDICAL PRODUCTS, INC
1485 Sarayashiki, Taraga-cho, Kasugai-city,
Aichi 486-0808 Japan
Tel: +81-568-81-7954
Fax: +81-568-81-7785

EC **REP**

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

TOKAI MEDICAL PRODUCTS

2019/08 (Rev.6)

ENGLISH TOKAI MICROCATHETER / INSTRUCTIONS FOR USE

DESCRIPTION

The Tokai Microcatheter "Carnelian" and "Carnelian Support" are catheters with a single lumen capable of accommodating guidewires which are smaller than the indicated size. A hub at the proximal end of the catheter incorporates a strain relief and a luer-lock adapter to facilitate the attachment of the accessories. A radiopaque marker at the distal tip facilitates fluoroscopic visualization. The outer diameter of the catheter is coated with a hydrophilic surface to enhance deliverability to the target vasculature.

INTENDED USE

The Tokai Microcatheter is intended to be used for super selective infusion of diagnostic, embolic, or therapeutic agents into peripheral vasculature (excluding central circulatory system and neurovasculature), and exchange/support of guidewires.

WARNINGS

- Contents are supplied sterile
- Do not use if pouch is opened or damaged. Inspect the package for any fault or damage to ensure the sterility of the microcatheter.
- Do not use the microcatheter in case of any surface irregularities, bends, or kinks. Any damage to the microcatheter may affect its performance and change its characteristics.
- Use the catheter before the Use By date indicated on the package.
- Single-use only
Do not re-use, re-process, or re-sterilize. Reuse, re-processing, or re-sterilizing may compromise the structural integrity and quality of the device resulting in patient injury, illness, or even death. Reuse, re-processing, or re-sterilizing may also create a risk of contamination of the catheter, cause infection or cross infection, or may be the source of an infectious disease being transmitted from one patient to another. Contamination of the catheter may lead to injury, illness, or death.
- Exercise care in handling of the catheter during a procedure. If any resistance is met, discontinue the procedure and identify the cause of the resistance by fluoroscopy. Movement of the microcatheter or guidewire against resistance may result in separation of the microcatheter or guidewire tip, damage to the microcatheter or guidewire tip, or vessel perforation.
- The infusion pressure with this catheter should not exceed 1000 psi. Infusion pressure in excess of the maximum pressure may result in catheter rupture or device failure.
Furthermore, Infusion pressure in excess of the maximum pressure may result in patients or users injury. If the flow of material through the catheter becomes restricted, do not attempt to clear the catheter lumen by infusion. Identify the cause of the blockage and take corrective action, or replace the catheter with a new one before resuming infusion.
- When infusing a drug or contrast medium, the catheter tip should be monitored under fluoroscopy. If no flow is observed, discontinue the infusion of materials, and replace the catheter with a new one. When the catheter lumen is obstructed, it may cause catheter expansion, rupture, or cracking resulting in leakage of materials even if the infusion pressure is lower than the maximum pressure.
- Comply with the instruction manuals for medical devices and drugs used together with the catheter.
- Before using the catheter, flush the lumen and outer surface with sterile heparinized saline to activate the hydrophilic coating. It is recommended that flush the lumen with sterile heparinized saline as needed to maintain lubrication during the procedure. If the catheter is used while the outer surface is poorly lubricated, the hydrophilic polymer coating will not function properly to maintain lubrication and the coating may come off during use.
- Overtightening of a hemostasis valve may cause damage to the catheter.

PRECAUTIONS

- If this catheter is inserted in a catheter with a stopcock, operating the stopcock may cause catheter rupture or device failure.
- Immersing this catheter in or wiping this catheter with chemicals containing organic solvents such as alcohol may cause damage to the catheter, or may compromise lubrication.
- Using the attached syringe for injecting oil-based contrast medium, fat emulsion, or preparations containing fat emulsion may cause damage to parts made of polycarbonate.
- When the catheter is used to replace or support the guidewire, do not rotate the tip of the catheter to prevent damage or disconnection.
- When using coiled embolization materials, select proper ones by checking the sizes of this product and aids in order to ensure that they fit together correctly.

<Precautions for Use>

- 1) Thoroughly flush the catheter with sterile heparinized saline solution prior to use.
- 2) Use the catheter immediately after the pouch is opened. Discard the product after use.
- 3) In case of a shaped product, be careful not to pinch the shaft with your fingers or the like strongly when pulling out the mandrel from the product. Also, if it is hard to pull out the mandrel, do not try to pull it out forcefully but gently little by little.
- 4) Handle with care to prevent accidental damage, bending, or kinking of the catheter when it is inserted into the guide catheter.
- 5) Exercise care in handling of the catheter during a procedure. If any resistance is met, discontinue the procedure and identify the cause of the resistance by fluoroscopy. Movement of the microcatheter or guidewire against resistance may result in separation of the microcatheter or guidewire tip, damage to the microcatheter or guidewire tip, or vessel perforation

COMPLICATIONS

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Possible complications include, but are not limited to, the following:

- Infection
- Perforation of vessel or arterial wall
- Angiospasm
- Hepatic dysfunction
- Non-target embolism
- Ischemia
- Hematoma
- Damage to vessel wall
- Artery dissection
- Embolization
- Death
- Inflammation
- Inflammatory contraction
- Allergic reaction
- Hepatic artery spasm
- Aneurysm
- Sepsis
- Abscess

INSTRUCTIONS FOR USE

- 1) Remove the microcatheter with hoop from the pouch.
- 2) Flush the hoop with heparinized saline to activate the hydrophilic coating. The luer fitting attached to the hoop may facilitate the flushing of the hoop.
- 3) Remove the microcatheter from the hoop, and confirm that the surface is well lubricated. If the microcatheter is difficult to remove, repeat injection or place in heparinized saline bath.
- 4) After removing the microcatheter from the hoop, thoroughly flush the lumen with sterile heparinized saline by attaching saline-filled syringe to the catheter hub. Use the attached syringe only for the purpose of cleaning the catheter lumen; do not use it for infusing drug or contrast medium into the body. Do not overtighten the syringe when connecting it to the hub.
- 5) Carefully Insert and advance an appropriate guidewire into the microcatheter.
- 6) Introduce the microcatheter wire assembly through a hemostatic valve. Carefully tighten the valve around the microcatheter to prevent backflow, allowing some movement through the valve by the microcatheter.
- 7) Advance the guidewire and the microcatheter to a selected vascular site by alternately advancing the guidewire and then tracking the microcatheter over the guidewire under fluoroscopy.
- 8) To infuse, completely remove the guidewire from the microcatheter. Connect a syringe with infusate or introduce embolic materials to the microcatheter hub and infuse as required.

CAUTION

This catheter should be used only by physicians thoroughly trained and experienced in percutaneous, intravascular techniques and procedures.

CONTRAINDICATIONS

Unknown

Tokai Medical Products, Inc. assumes no liability with respect to devices reused, re-processed, or re-sterilized. Tokai Medical Products, Inc. makes no warranties, expressed or implied, including, but not limited to, merchantability and fitness to the intended use with respect to such device.

Any serious or life-threatening adverse events or deaths associated with use of this microcatheter should be reported to the device manufacturer and Competent Authority of the country where it occurred.

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