# Fortress

---

Introducer Sheath System

Einführbesteck System

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Sistema de Funda de Introducator

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Sistem vodila

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**STERILE EO**

**Caution: Federal (USA) law  
restricts this device to sale  
by or on the order of a  
physician.**

**BIOTRONIK**

excellence for life

---

<b>English</b>			
Instructions for Use	Page		4
<b>Deutsch</b>			
Gebrauchsanweisung	Seite		6
<b>Français</b>			
Mode d'utilisation	Page		8
<b>Italiano</b>			
Istruzioni per l'uso	Pagina		10
<b>Español</b>			
Instrucciones de uso	Página		12
<b>Български</b>			
Инструкции за употреба	стр.		14
<b>Hrvatski</b>			
Upute za uporabu	stranica		16
<b>Čeština</b>			
Návod k použití	Strana		18
<b>Dansk</b>			
Brugsanvisning	Side		20
<b>Nederlands</b>			
Gebruiksaanwijzing	Pagina		22
<b>Suomi</b>			
Käyttöohjeet	Sivu		24
<b>Ελληνικά</b>			
Οδηγίες Χρήσης	σελίδα		26
<b>Magyar</b>			
Használati utasítás	Oldal		28
<b>Latviski</b>			
Lietošanas instrukcija	lpp.		30
<b>Lietuviškai</b>			
Naudojimo instrukcijos	p.		32
<b>Norsk</b>			
Bruksanvisning	Side		34
<b>Polski</b>			
Instrukcja użytkowania	Strana		36
<b>Português</b>			
Instruções de utilização	Página		38
<b>Română</b>			
Instrucțiuni de utilizare	pagina		40
<b>Русский</b>			
Инструкция по применению	Страница		42
<b>Slovensky</b>			
Návod na použitie	strana		44
<b>Slovenščina</b>			
Navodila za uporabo	Stran		46
<b>Svenska</b>			
Bruksanvisning	Sidan		48
<b>Türkçe</b>			
Kullanım Kılavuzu	Sayfa		50

---

**Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.**

### Product description

The Fortress Introducer Sheath System consists of an introducer sheath with haemostasis valve and side port, as well as dilator with a tapered tip and luer lock at the proximal end. The main introducer sheath tubing is connected at the proximal end to a haemostasis valve with side port tubing that is connected to a plastic 3-way stopcock valve. The side port is used for flushing the introducer sheath. The introducer sheath is introduced into the vascular system with the aid of the dilator. The haemostasis valve at the proximal end of the introducer sheath conforms and seals around guide wires and catheters to reduce blood leakage from the introducer sheath. A radiopaque marker helps identify the distal end of the introducer sheath. The 5F (1.95 mm), 6F (2.28 mm), 7F (2.55 mm) and 8F (2.88 mm) Introducer Sheath has a hydrophobic silicone coating on the outer surface of the distal 30 cm portion.

The Fortress Introducer Sheath System consists of the following components:

- One (1) introducer Sheath with haemostasis valve
- Two (2) dilators; 4F (1.62 mm) 45 cm and 100 cm; 5F (1.95 mm) and 6F (2.28 mm) 100 cm
- One (1) dilator; 4F (1.62 mm) 65 cm and 90 cm; 5F (1.95 mm), 6F (2.28 mm), 7F (2.55 mm) and 8F (2.88 mm) 45 cm, 65 cm and 90 cm

### Indications

The Fortress Introducer Sheath System is intended to provide access and to facilitate percutaneous introduction of guide wires, catheters and other devices into the femoral, popliteal and infrapopliteal arteries while maintaining haemostasis during diagnostic and interventional procedures.

### Contraindications

- The Fortress Introducer Sheath System is not designed, sold or intended for any use other than as indicated.
- It is the responsibility of the physician to determine whether any physical impairment of the patient would contraindicate the use of this device.
- Do not attempt introducer sheath insertion or use if the following conditions apply:
  - Infection or lesion near the insertion site
  - An adjacent vascular access device that could compromise catheter care and cleanliness
  - Anatomical irregularities in the patients extremities that could interfere with proper placement of the catheter
  - Prior trauma to vessels that could interfere with insertion of the catheter

### Warnings

- Read the instructions prior to use.
- Do not use the device if the package or contents appear damaged in any way.
- This device is intended for single use only. Do not reuse, resterilize or autoclave. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. Contract Medical International GmbH (CMI) will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- Do not attempt to advance or withdraw the introducer, guide wire, catheter, or other interventional device if resistance is felt. Use fluoroscopy to determine the cause. If the cause cannot be determined and corrected, discontinue the procedure and withdraw the introducer sheath. Continued advancement or retraction against resistance may result in serious injury, and/or breakage of the guide wire, introducer sheath, catheter or interventional medical device.
- Do not use the device with a power injector.
- Do not alter this device. Alterations may impair function. Be aware that attempting to alter the shape of a coated device by twisting or similar methods may compromise the coating integrity and the damage to the coating could result in coating material remaining in the vasculature, which may require additional intervention.
- Verify introducer sheath, device, catheter and accessory components compatibility prior to use.
- Protection tube from sheath and dilator must be removed prior to use.
- Do not attempt to use a guide wire with a maximum diameter greater than specified on package label. Device damage or breakage may occur.
- Do not attempt to insert a catheter or other interventional device with a diameter larger than the introducer sheath size indicated. Device damage or breakage may occur.
- This device should only be used by or under the direction of physicians thoroughly trained in the technique of catheter delivery systems.
- Observe sterile technique at all times when handling and inserting or removing the catheter.
- Adequate heparinization must be maintained during the procedure. The risk of systemic anticoagulation must be weighed against the benefits of extracorporeal circulation.
- The medical techniques and procedures described in these instructions are presented as an example only, and do not represent ALL medically acceptable protocols. They are not intended as a substitute for the physician's experience and judgment in treating any specific patient.
- Device contains metal components: do not expose to MRI.

### Precautions

- Store in a cool dry place.
- Follow manufacturer's instructions for proper device storage as improper storage can impact the integrity of coating.
- Do not expose to organic solvents.
- Examine packaging and device before use. Do not use if either the packaging or device is damaged or if sterile barrier has been compromised.
- The device is supplied sterile. The package will serve as an effective barrier until the "use by" (expiration) date printed on the box.
- Do not attempt introducer sheath advancement or withdrawal without guide wire and dilator secured in place. Severe vascular damage and/or injury may occur.
- Insertion into and removal from artery may cause excessive bleeding and/or other complications.
- Individual patient anatomy and physician technique may require procedural variations.

### Potential Complications

Potential complications include, but are not limited to:

- Air embolism
- Bleeding
- Vessel spasm
- Haematoma
- Infection/ Sepsis
- Thromboembolism
- Thrombophlebitis
- Thrombosis
- Perforation or laceration of the vessel wall
- Intimal tear
- Risks normally associated with percutaneous diagnostic and/or interventional procedures

Some of the complications listed above could result in death or serious injury to the patient.

## Directions for use

### INSERTION:

1. Prepare the insertion site using aseptic technique and local anaesthesia as necessary.
2. Remove the Fortress introducer sheath and dilator from the packaging using aseptic technique and examine for possible damage and defects. Do not use in case of visible damage.
3. If the product has a detachable haemostasis valve, ensure this is securely tightened to the hub of the sheath.
4. Flush dilator, introducer sheath and side port with saline, suitable isotonic solution or heparinized intravenous fluid. After flushing side port, turn stopcock to off position to keep flushing solution in side port and prevent bleeding after insertion into the vessel.
5. Insert dilator tip through the haemostasis valve of the Fortress Introducer Sheath System completely into the introducer sheath until the dilator hub comes in contact with the haemostasis valve to ensure that the tapered portion of the dilator is beyond the end of the introducer sheath. When advancing the introducer sheath, ensure that the dilator remains fully inserted into the introducer sheath to avoid damage to the vessel.

**Note:** When using product supplied with the pre-shaped dilator, use extra care when inserting this pre-shaped dilator into the introducer sheath as the pre-shaped dilator increases the friction.

6. Follow standard recommended practice for vessel puncture or incision. While holding the access needle, place flexible or J-end of the guide wire through access needle into the vessel.

**Note:** Refer to product labelling for the guide wire size compatible with system components.

**Warning:** Do not advance the guide wire if resistance is met. Determine cause of resistance before proceeding.

7. Hold guide wire in place and remove access needle. Hold pressure at the site until the introducer sheath/dilator assembly is placed.
8. Insert the introducer sheath/dilator assembly over the guide wire into the vessel. Under fluoroscopic guidance, carefully advance the introducer sheath/dilator assembly over the guide wire until it is at the desired location. Do not allow dilator to back out of the introducer sheath while advancing. Stop advancement if resistance is met and investigate cause before proceeding.
9. After the introducer sheath/dilator assembly has been positioned to the desired location, hold the introducer sheath steady and maintain guide wire position while withdrawing the dilator from the introducer sheath over the guide wire until it is completely removed.
10. Aspirate from the side port extension to remove any potential air. After aspiration, flush the side port with a suitable solution. Stopcock should be turned off to maintain fluid in side port.
11. While maintaining the position of the guide wire, advance catheter or other interventional device over the guide wire into the introducer sheath.

**Note:** Hold the introducer sheath in place when inserting, positioning, or removing the devices.

**Note:** When the introducer sheath will remain in a vessel for an extended time period, consider using a continuous drip of heparinized intravenous fluid under pressure administered through the side port connection.

### REMOVAL:

1. Insert the dilator over the guide wire and fully into the introducer sheath.
2. The introducer sheath may be removed when clinically indicated. Compression on the vessel above the puncture site should be started as the introducer sheath is slowly removed.

**Note:** Collected fibrin at the tip of the introducer sheath may be aspirated via the side arm tubing prior to removal of the introducer sheath.

3. Discard the introducer sheath appropriately.
4. Upon removal of the introducer sheath, precautions should be taken to prevent bleeding, vessel damage, or other serious injury.

## Symbol legend

Explanation of symbols used on the package labels:

Single use only

Sterilized using ethylene oxide

See instructions for use before use

Do not re-sterilize

Reference number

Do not use if package is damaged

Lot number

Keep away from sunlight

Use by Date

Manufacturer

Temperature limitation

Date of manufacture

Keep Dry

## Warranty / Liability

The product and each component of its system (hereinafter „the product“) have been designed, manufactured, tested and packaged with all reasonable care. However, CMI has no control over the conditions under which the product is used and a disturbance of the intended function of the product may occur for various reasons. In this respect, the warning in this product publication/instructions for use are expressly to be considered as an integral part of this Disclaimer and provide more detailed information. For this reason, CMI disclaims expressly all warranties, expressed or implied regarding the product, including but not limited to, any warranty of merchantability or fitness for a particular purpose of the product. Product descriptions or user guidelines in publications do not constitute any expressed representation or any expressed warranty. CMI is not liable for any direct, incidental or consequential damages or medical expenses caused by any use, defect, failure or malfunction of the product whether the claim is based on contract, warranty, tort or otherwise.

This does not apply in the case of intention or in the case of gross negligence of legal representatives or executive staff of CMI. In commercial transactions relating to merchants, the liability is limited to the compensation of typical damages: compensation for any untypical or incidental damage is excluded. These limitations of liability and warranty are not intended to contravene any mandatory provision of law applicable in the respective country. If any clause of the Disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects CMI's legitimate interest in limiting its liability or warranty without infringing any mandatory provisions of applicable law. No person has any authority to bind CMI to any warranty or liability regarding the product.

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All rights reserved. Specifications are subject to modification,
BIOTRONIK AG