English	Peripheral Multifunctional Catheter	Instructions for Use	9	en
Deutsch	Peripherer Multifunktionskatheter	Gebrauchsanweisung	10	de
Français	Cathéter périphérique multifonctionnel	Mode d'emploi	12	fr
Italiano	Catetere periferico multifunzionale	Manuale tecnico di istruzioni	14	it
Español	Catéter multifuncional periférico	Instrucciones de uso	15	es
Български	Периферен мултифункционален катетър	Ръководство за употреба	17	bg
Hrvatski	Periferni višenamjenski kateter	Priručnik za uporabu	19	hr
Česky	Periferní multifunkční katétr	Návod k použití	20	cs
Dansk	Perifert multifunktionelt kateter	Brugsanvisning	22	da
Nederlands	Perifere Multifunctionele Katheter	Technische handleiding	23	nl
Eesti	Perifeerne multifunktsionaalne kateeter	Kasutusjuhend	25	et
Suomi	Perifeerinen monitoimikatetri	Käyttöohje	27	fi
Ελληνικά	Περιφερικός πολυλειτουργικός καθετήρας	Οδηγίες χρήσης	28	el
Magyar	Perifériás multifunkciós katéter	Használati utasítás	30	hu
Latviešu	Perifērs daudzfunkcionāls katetrs	Lietošanas norādījumi	31	lv
Lietuvių	Periferinis daugiafunkcinis kateteris	Naudojimo instrukcija	33	lt
Norsk	Perifert multifunksjonelt kateter	Bruksanvisning	35	no
Polski	Wielofunkcyjny cewnik obwodowy	Instrukcja obsługi	36	pl
Português	Cateter periférico multifuncional	Manual técnico	38	pt
Română	Cateter periferic multifuncțional	Instrucțiuni de utilizare	40	ro
Русский	Периферический многофункциональный катетер	Инструкция по применению	41	ru
Slovenčina	Periférny multifunkčný katéter	Návod na použitie	43	sk
Slovenščina	Periferni večfunkcijski kateter	Navodila za uporabo	45	sl
Svenska	Perifer multifunktionskateter	Bruksanvisning	46	sv
Türkçe	Periferik Çok Fonksiyonlu Kateter	Kullanım Talimatları	48	tr

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Description

The Oscar® Peripheral Multifunctional Catheter system (hereinafter Oscar or the device) is intended for dilatation of stenotic segments in peripheral vessels. The device is intended to guide and support a guidewire during access and crossing of lesions in the peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents. Oscar encompasses three components: Oscar Support Catheter (Figure 1) with integrated Lock Grip (Figure 1.2), Oscar Dilator (Figure 2) and Oscar PTA Balloon (Figure 3). Single Oscar PTA Balloon is also available separately, and can only be used together with the pre-supplied Oscar Support Catheter as indicated in the Compatibility Chart. Oscar Support Catheter and Oscar Dilator are used in tandem (Figure 4) to facilitate lesion access and lesion crossing by guide wire. The Oscar Dilator is removed from the Oscar Support Catheter after crossing a lesion with a compatible guide wire. Once the Oscar Dilator is removed, the Oscar PTA Balloon is inserted through the Oscar Support Catheter over a compatible guide wire for balloon dilatation (Figure 5).

1. Oscar Support Catheter with integrated Lock Grip

The Oscar Support Catheter is a single-lumen support catheter offered in 4F and 6F introducer sheath compatible configurations (see Compatibility Chart). The 4F configuration is available with a usable length of 108 cm (long). The 6F configuration is available with 60 cm (short) or 108 cm (long) usable length shaft. The Oscar Support Catheter includes a radiopaque marker (Figure 1.1) located on the distal tip. The shaft is hydrophobically coated over the distal 45 cm. Located at the proximal end of the Oscar Support Catheter is a Lock Grip (Figure 1.2) that has three parts. The Flush Port (Figure 1.3) is used with a regular three-way stopcock and intended for injection of fluids including sterile saline or contrast media. Next to the Flush Port is the Tip-Control (Figure 1.4) which has a mechanical, spring-loaded actuation mechanism, and the Seal+Lock Cap (Figure 1.5) which acts as a lockable hemostatic valve.

The Tip-Control is used to extend the Oscar Dilator to allow injection of fluids without the need to remove the dilator or guide wire.

The Seal+Lock Cap is closed to seal and secure the respective positions of either Oscar Dilator, Oscar PTA Balloon or a guide wire only during use. When the Seal+Lock Cap is in the open position, the Lock Grip permits the manual insertion, advancement, positioning and removal of the Oscar Dilator or Oscar PTA Balloon and a guide wire.

2. Oscar Dilator

The Oscar Dilator is used together with the Oscar Support Catheter. The radiopaque tip (Figure 2.1) of the Oscar Dilator is aligned with the distal end of the Oscar Support Catheter in the so called "0" [Zero] position (Figure 7.2.1) and tapers down for atraumatic access. The Oscar Dilator can be extended up to 19.1 cm from the Oscar Support Catheter and is compatible with guide wires according to the label.

3. Oscar PTA Balloon

The Oscar PTA Balloon is an over-the-wire (OTW) percutaneous transluminal angioplasty (PTA) dilatation catheter, designed for variable-length balloon inflations (Figure 3.1) with evenly spaced radiopaque markers (Figure 3.2) to facilitate visualization and positioning when used in combination with the Oscar Support Catheter. The dilatation balloon is designed to be inflated to a known diameter at a specific inflation pressure consistent with the Compliance Chart.

The Oscar PTA Balloon includes a soft tapered tip (Figure 3.3) to facilitate advancement of the catheter. The inflation/deflation port (Figure 3.4) serves for connecting an inflation device to inflate/deflate the balloon. The guide wire port (Figure 3.5) enables flushing of the guide wire lumen. The inner surface of the guide wire lumen and the outer surface of balloon tip to distal balloon cone are hydrophobically coated.

When the balloon is positioned across the lesion, the balloon length is adjusted by distally exposing a desired balloon length from the distal tip of the Oscar Support Catheter. The exposed length is determined as the distance between the most distal balloon marker and the radiopaque marker of the Oscar Support Catheter. Once the length of the balloon has been adjusted, the Oscar PTA Balloon position is locked for lesion dilatation using the Seal+Lock Cap. The balloon can be inflated and deflated until no contrast medium is visible in the balloon and withdrawn into the Oscar Support Catheter aligning distal balloon marker with the Oscar Support Catheter radiopaque marker, for either subsequent lesion dilatation or PTA balloon is removed entirely alone or together with the Oscar Support Catheter.

The Oscar PTA Balloon is compatible with guide wires and the Oscar Support Catheter according to the label.

How Supplied

Sterile. Non-pyrogenic. The device is sterilized with ethylene oxide.

Contents

The device is supplied in two trays within a sealed, peel-open pouch. The device is supplied with one Instruction for Use manual. The trays contain the following components:

- Tray Support: one Oscar Support Catheter and one Oscar Dilator
 - Tray PTA: one Oscar PTA Balloon and one Compliance Card
- Each Oscar PTA Balloon is packaged with a protective tubing that is positioned over the balloon.

Single Oscar PTA Balloons are also separately available and can only be used together with the pre-supplied Oscar Support Catheter (please read the Compatibility Chart).

Storage

Store protected from sunlight and keep dry. Store between 15°C [59°F] and 25°C [77°F]. Short-term storage temperature excursions between 10°C [50°F] and 40°C [104°F] are allowed up to 48 hours.

Indications

Oscar is indicated for percutaneous transluminal interventions in the peripheral vasculature to provide support during access into and to dilate stenoses in femoral, popliteal and infrapopliteal arteries. The product is also intended for injection of radiopaque contrast media for the purpose of angiography.

Contraindications

All general contraindications for percutaneous transluminal angioplasty (PTA) are contraindications for this device. Contraindications for this device and peripheral dilatation catheters in general are:

- Lesions that cannot be reached or treated with the system
- Uncorrected bleeding disorders
- Sepsis

Furthermore, all general PTA and procedure-related contraindications as described in the national and international guidelines of the respective medical associations apply.

Warnings

- Oscar has one intended use which can only be fulfilled by using all components together, Oscar Support Catheter with Oscar Dilator or Oscar PTA Balloon, and the components are not compatible with other devices or must not be used for another purpose. DO NOT use regular (non-Oscar) PTA balloons in combination with the Oscar Support Catheter.
- Before use, visually inspect the packaging. If the packaging is damaged or opened, or if any information is obscured, DO NOT use the device.
- DO NOT use the device in the coronary, cervical and intra-cranial arteries as this is prohibited.
- Open products must be disposed after the procedure due to compromised sterility.
- The device is designed and intended for single use only. DO NOT resterilize and/or reuse. Reuse of single-use devices creates a potential risk of patient or user infections. Contamination of the device may lead to injury, illness or death of the patient. Cleaning, disinfection and sterilization may compromise essential material and design characteristics leading to device failure. BIOTRONIK will not be responsible for any direct, incidental or consequential damages resulting from resterilization or reuse.
- The inflated diameter of the balloon should never exceed the original diameter of the vessel proximal and distal to the lesion.
- DO NOT exceed the rated burst pressure (RBP) indicated on the Compliance Chart. The use of a pressure-monitoring device is mandatory to prevent over-pressurization.
- When the system is introduced into the vascular system, it should be manipulated under high-quality fluoroscopy.
- Before advancement of the Oscar Support Catheter, the pre-inserted Oscar Dilator or Oscar PTA Balloon are required to be positioned in the "0" position, with their respective tips aligned to the Oscar Support Catheter tip and locked in place using the Seal+Lock Cap.
- DO NOT use the device after the use by date indicated on the label.
- DO NOT use air or any gaseous medium to inflate the balloon. Use only an appropriate balloon inflation medium (e.g. 50:50 mixture by volume of contrast medium and sterile saline).
- DO NOT use the device if a vacuum cannot be held as this indicates a leak in the system.
- Ensure the Seal+Lock Cap is in the closed position before inflating the Oscar PTA Balloon or performing sterile saline/contrast medium injection with the Oscar Dilator.
- DO NOT expose the device to organic solvents e.g. alcohol.
- DO NOT use a power injection system when the Dilator is inserted into the Support Catheter.

Precautions

General Precautions

- The device should only be used in medical facilities by physicians thoroughly trained and experienced in performing vascular interventions, (including cases of life-threatening complications).
- Prior to the procedure, the device should be visually examined for any damage to verify its functionality and to ensure that its size is suitable for the specific procedure for which it is to be used.
- Use only guide wires compatible with the 4F [0.014"] and 6F [0.018"] configurations.
- In order to prevent and reduce clotting, flush or rinse all products entering the vascular system with sterile saline or a similar solution prior to use. The use of systemic heparinization during the procedure is recommended.
- Non-clinical testing with Cordis AVANTI®, Cordis BRITE TIP™, Cook's Flexor® Check-Flo®, Teleflex' Super Arrow-Flex® and BIOTRONIK's Fortress introducer sheaths has demonstrated that the device is compatible with the indicated minimum introducer sheath sizes. If the device is used with other coiled, braided or non-braided sheaths, a larger French size than indicated on the label may be necessary to reduce friction.
- Exercise care during handling to reduce the possibility of accidental breakage, bending or kinking of the device.
- If a strong resistance is felt during manipulation, stop the procedure and determine the cause before proceeding.
- DO NOT torque Oscar Support Catheter or Oscar Dilator excessively if it is already bent, in order to avoid separation or breakage of the product.
- Non-clinical testing has shown Oscar PTA Balloon re-folding and length adjustment capability for three (3) inflations up to RBP.
- Non-clinical testing has shown appropriate functionality of the Oscar Support Catheter for six (6) inflations of Oscar PTA Balloon up to RBP.
- Constituents of the device contain the following substances defined as CMR 1B in a concentration above 0.1% w/w:
 - Cobalt: CAS No. 7440-48-4; EC No. 231-158-0
 - 6,6'-di-tert-butyl-2,2'-methylene-di-p-cresol: CAS No. 119-47-1; EC No. 204-327-1
- Current scientific evidence supports that medical devices manufactured from cobalt alloys or stainless steel alloys containing cobalt DO NOT cause an increased risk of cancer or adverse reproductive effects. Based on literature-based toxicological evaluation, exposure to the amount of CAS 119-47-1 contained in the device DOES NOT cause an increased risk of adverse reproductive effects.

Potential Adverse Events/Complications

Possible complications include, but are not limited to:

- Injury to the vessel wall, intimal tear
- Arteriovenous fistula
- Embolization of air, thrombotic or atherosclerotic material
- Pseudoaneurysm formation
- Vessel spasm
- Dissection
- Restenosis of the dilated vessel

- Total occlusion of the vessel
- Infection
- Hemorrhage or hematoma
- Embolism
- Thrombosis
- Abrupt closure
- Allergic reactions to contrast media, antiplatelets, anticoagulants
- Catheter events: failure to reach or cross the lesion, inflation difficulties, rupture or pinhole of the balloon, deflation difficulties, withdrawal difficulties, embolization of catheter material
- Emergency surgery to correct vascular complications
- Tissue necrosis and limb loss
- Perforation and rupture of the vessel
- Nephropathy
- Death.

Furthermore, all procedure-related events as described in the national and international guidelines of the respective medical associations apply.

Directions for Use

Oscar Dilator and Oscar PTA Balloon can only be used in combination with Oscar Support Catheter.

1. Oscar Support Catheter and Oscar Dilator

Device Preparation

1. Check and ensure that the outer packaging is not damaged or opened prior to use. Open the outer packaging and remove the pouch.
2. Inspect the pouch and ensure that there is no damage to the sterile barrier before proceeding. If loss of sterility is suspected, DO NOT use the device. Using an aseptic technique, open the pouch and remove the sterile content.
3. Remove the Oscar Support Catheter/Oscar Dilator combination from Tray Support.

Flushing and Air Removal

Note: DO NOT remove the Dilator from the support catheter.

4. Connect a three-way stopcock to the Flush Port of the Oscar Support Catheter.
5. Connect a 10 ml or 20 ml syringe containing sterile saline to the three-way stopcock on the Flush Port of the Oscar Support Catheter.
6. Confirm that the Seal+Lock Cap of the Oscar Support Catheter is in the open position (Figure 7.1.1) by moving the Oscar Dilator.
7. Flush the Oscar Support Catheter until fluid visibly exits from the proximal end of the Seal+Lock Cap.
8. Turn the Seal+Lock Cap into the closed position (Figure 7.1.2), confirmed by no movement of the Oscar Dilator.
9. Flush the Oscar Support Catheter until fluid visibly exits from the Support Catheter distal end. After flushing, turn the three-way stopcock to the off position.
10. Open the Seal+Lock Cap to align the Oscar Dilator tip with the distal tip of the Oscar Support Catheter to position "0". Close the Seal+Lock Cap of the Oscar Support Catheter to secure the position of the Oscar Dilator and confirm visually that the tips are aligned.
11. Attach a syringe to the Oscar Dilator hub and flush the guide wire lumen until fluid visibly exits the distal end of the Oscar Dilator and air has been removed completely.

Note: Closing the Seal+Lock Cap is only permissible on the metal shaft of either Oscar PTA Balloon or Oscar Dilator.

Insertion/Use

12. Insert Oscar Support Catheter and Oscar Dilator as a single unit over a compatible guide wire and advance them to the target lesion under fluoroscopy.
13. To adjust the level of guide wire support, open the Seal+Lock Cap and advance or retrieve the Oscar Dilator.
14. After successful crossing, gently withdraw the Oscar Dilator from the Oscar Support Catheter by opening the Seal+Lock Cap.
15. Close the Seal+Lock Cap after Dilator removal to avoid bleeding.

Saline/Contrast Injection

Oscar Support Catheter is designed to enable sterile saline and/or contrast media injection.

Saline and/or Contrast Injection through Oscar Support Catheter without removing Oscar Dilator

TIP-CONTROL MODE:

1. Maintain the Seal+Lock Cap of the Oscar Support Catheter in the closed position.
2. Mechanically actuate the Oscar Dilator tip distally by performing the following steps (Figure 7.2):
 - Turn the Tip-Control of the Lock Grip (Figure 7.2.1) to its counter clockwise position to unlock the mechanical actuation mode (Figure 7.2.2).
 - Move the Tip-Control to its most distal position, and turn the Tip-Control clockwise to lock the position (Figure 7.2.3).
3. Perform contrast medium injection via the Flush Port of the Oscar Support Catheter.
4. Retract the Oscar Dilator tip back in access position (Oscar Dilator tip aligned with Oscar Support Catheter tip) by performing the following steps:
 - Turn the Tip-Control counter clockwise to unlock the position,
 - Release the Tip-Control to its most proximal position, and
 - Turn the Tip-Control clockwise to lock the position.

MANUAL MODE:

1. Open the Seal+Lock Cap of the Oscar Support Catheter.
2. Manually extend the Oscar Dilator distally from the Oscar Support Catheter by a minimal distance of 2 mm under fluoroscopic guidance.
3. Close the Seal+Lock Cap to secure the Oscar Dilator in position.
4. Perform contrast medium injection via the Flush Port of the Oscar Support Catheter.
5. Open the Seal+Lock Cap; realign Oscar Dilator tip to Oscar Support Catheter tip; and close the Seal+Lock Cap for subsequent advancement of the Oscar Support Catheter.

Note: Before injecting contrast media, connect the syringe to the three-way stopcock on the Flush Port of the Oscar Support Catheter.

Note: Ensure that the three-way stopcock is in the open position in order to inject the contrast media.

Note: After flushing, turn the three-way stopcock to the off position.

Saline and/or contrast injection through Oscar Support Catheter without Oscar Dilator

Caution: The use of a power injection system can lead to device failure. Limit the injection pressure to 300 psi (refer to Flow Rates table for contrast media flow rates).

1. Open the Seal+Lock Cap of the Oscar Support Catheter.
2. Remove the Oscar Dilator entirely from the Oscar Support Catheter while maintaining position of the guide wire.
3. Close the Seal+Lock Cap.
4. Perform contrast medium injection via the Flush Port of the Oscar Support Catheter.
5. Open the Seal+Lock Cap.
6. Reinsert the Oscar Dilator over the guide wire.
7. Realign the Oscar Dilator tip to the Oscar Support Catheter tip into position "0" under fluoroscopic guidance.
8. Close the Seal+Lock Cap.

Note: DO NOT advance the Oscar Support Catheter without the Oscar Dilator.

Note: Before injecting contrast media, connect the syringe or injection system to the three-way stopcock on the Flush Port of the Oscar Support Catheter.

Note: Ensure that the three-way stopcock is in the open position in order to inject the contrast media.

Note: After flushing, turn the three-way stopcock to the off position.

2. Oscar PTA Balloon

Device Preparation

1. Remove the Oscar PTA Balloon from Tray PTA, then remove protective tubing from the Oscar PTA Balloon.
2. Connect a 10 ml or 20 ml syringe containing sterile saline to the guide wire port.
3. Flush the guide wire lumen.
4. Remove the syringe.

Purge air from catheter inflation lumen

5. Fill an inflation device with contrast medium solution.
6. Remove air from the inflation device according to the manufacturer's recommendations and instructions.
7. Connect a three-way stopcock to the inflation/deflation Luer port of the Oscar PTA Balloon and attach the inflation device to the three-way stopcock.
8. Open the three-way stopcock to create an open fluid path between the Oscar PTA Balloon and the inflation device.
9. Draw vacuum for 30 seconds to remove air.
10. Close the stopcock; purge the inflation device of all air.
11. Repeat steps 8 to 10 to purge air.
12. Open the stopcock; pull vacuum and hold.

Insertion

1. Verify that the Oscar Support Catheter is in place and positioned proximal to the lesion.
2. Activate the Tip-Control by turning it counterclockwise, then push it distally and fixate it by turning the Tip-Control clockwise.
3. Connect a 10 ml or 20 ml syringe containing sterile saline to the three-way stopcock on the Flush Port of the Oscar Support Catheter and flush.
4. Open the Seal+Lock Cap of the Oscar Support Catheter.
5. Insert the Oscar PTA Balloon over the guide wire into the Oscar Support Catheter until the distal balloon marker is aligned with the radiopaque marker of the Oscar Support Catheter.

The balloon with dimensions 2 to 6 mm in diameter has 3 markers starting distally, each 60 mm apart. The 7 mm balloon has 2 markers starting distally, 50 mm apart.

6. Position the Oscar PTA Balloon across the lesion by exposing the desired balloon length. The exposed length is determined as the distance (Figure 6) between the radiopaque marker of the Oscar Support Catheter and the most distal balloon marker. The exposed balloon length should be at least 20 mm. Verify fluoroscopically that the balloon length exposed is adequate for the length of the lesion to be treated.

Note: Alternatively, if the Oscar Support Catheter has crossed the lesion, pull the Oscar Support Catheter back until the desired balloon length has been exposed.

Note: If the Oscar PTA Balloon is intended to be used at maximum length of 180 mm (balloons from 2 to 6 mm in diameter) or 100 mm (balloon with 7 mm diameter), leave about 5 mm of the metal shaft of the Oscar PTA Balloon visible.

Inflation/Deflation

8. Inflate the balloon to achieve the desired dilatation using standard PTA techniques.
9. Deflate the balloon by pulling vacuum with the inflation device until contrast medium is no longer visible in the balloon.

Caution: Verify under fluoroscopy if the contrast medium has been removed from the balloon. Longer and larger balloons require typically longer deflation time.

10. Open the Seal+Lock Cap.

11. Pull the Oscar Support Catheter about 5 mm back while maintaining vacuum. Then withdraw the Oscar PTA Balloon into the Oscar Support Catheter until the distal radiopaque marker of Oscar PTA Balloon is aligned with the Oscar Support Catheter marker.

12. For dilatation with an adjusted balloon length open the Seal+Lock Cap of the Oscar Support Catheter and repeat steps 6-11.

Note: During adjustments of the Oscar PTA Balloon ensure that the metal shaft is always inside the Oscar Support Catheter.

Note: Another Oscar PTA Balloon compatible (please read the Compatibility Chart) with the Oscar Support Catheter may be used if needed.

Note: Before re-adjusting the Oscar PTA Balloon length, the balloon is required to be fully deflated with no contrast media visible in the balloon and the Seal+Lock Cap of the Oscar Support Catheter unlocked.

Removal

13. Following completion of the PTA, deflate the balloon, hold vacuum and withdraw the Oscar PTA Balloon into the Oscar Support Catheter until the distal balloon marker is aligned with the radiopaque marker of the Oscar Support Catheter. In case another Oscar PTA Balloon is planned to be used, remove the used Oscar PTA Balloon while maintaining guide wire position and follow the steps described in the Oscar PTA Balloon Section. Otherwise lock the Oscar PTA Balloon and remove both components together while maintaining guide wire position.

Note: If the Oscar PTA Balloon cannot be completely withdrawn into the Oscar Support Catheter, then withdraw the Oscar PTA Balloon together with the Oscar Support Catheter.

14. After use, the device and packaging may be contaminated with infectious substances (e.g. blood). Therefore, dispose of the device and packaging in accordance with hospital, administrative and/or local government policy.

Warranty/Liability

This product and each of its components (hereinafter product) were designed, manufactured, tested and packaged with all reasonable care. However, since BIOTRONIK does not have control over the conditions under which the product is used, the contents of this Instructions For Use (IFU) are to be considered as an integral part of this disclaimer for cases when a disturbance of the product's intended function may occur for various reasons.

BIOTRONIK does not guarantee that the following events will not occur:

- Product malfunctions or failures
 - Patient's immune response to the product
 - Medical complications during the use of the product or as a consequence of the product being in contact with the patient's body
- BIOTRONIK does not assume any liability for:
- The use of the product that is not in accordance with the stated intended use/indication, contraindications, warnings, precautions and the directions for use of this IFU
 - Modification to the original product
 - Events which could not have been foreseen at the time of product delivery using the available levels of science and technology
 - Events originating from other BIOTRONIK products or products not from BIOTRONIK
 - Force majeure events including, without being limited to, natural disasters.

The above provisions shall be without prejudice to any disclaimer and/or limitation of liability agreed separately with the customer to the extent permitted by applicable laws.

Intellectual Property

The following list includes trademarks or registered trademarks of the BIOTRONIK Group of Companies filed in the EU, the United States, and possibly in other countries: BIOTRONIK, Oscar and Fortress. All other trademarks are the property of their respective owners.

Clinical Benefits

The use of the device for the described indications is expected to result in lumen enlargement leading to an improved blood flow of the treated vessel segment and subsequent perfusion of the lower extremities, thereby allowing relief from symptoms and improving quality of life, as well as providing a safe and effective intervention.

Reporting of Serious Incidents

Report any serious incident that has occurred with the device to the manufacturer and to the competent authority of the country in which your hospital is established. In case of a device malfunction, stop using the device immediately and return it to the manufacturer.

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All rights reserved. Specifications subject to
BIOTRONIK AG
© 2023 BIOTRONIK AG