



12-month results

## Conclusions

- Application of a paclitaxel coated balloon using Butyryl-tri-hexyl citrate (BTHC) as excipient is feasible and safe in a mixed population of patients with predominantly type I BMS or DES-ISR lesion
- A short exposure of the vessel wall to paclitaxel results in very low Late Lumen Loss (LLL), revascularization and MACE rates
- Pantera Lux application is a valuable treatment option for In-Stent Restenosis (ISR) in either BMS or DES patients

## Study design

Prospective, multi-center, non-randomized, European clinical trial of the Pantera Lux Paclitaxel Releasing Balloon

## Endpoints

### Primary endpoint

- LLL at 6 months

### Secondary endpoint

- MACE\* at 6 and 12 months


## Inclusion/exclusion criteria

### Major inclusion criteria

- Patients with a single restenotic lesion
- Target reference vessel diameter 2 – 4 mm
- Target lesion length 8 – 28 mm
- Target lesion stenosis  $\geq 50\%$  –  $< 100\%$

### Major exclusion criteria

- Myocardial Infarction (MI)
- Additional coronary lesions in the same vessel which requires treatment
- Totally occluded coronary artery



## Patient characteristics

Demographics	
Age (years)**	66.1 ± 9.4
Male gender	77.8%
Medical history/risk factors	
Hyperlipidemia	87.7%
Hypertension	87.7%
Prior MI	63.0%
Diabetes	27.2%
Renal disease	13.6%


## Lesion characteristics

Mehran classification	
Focal (Type I)	71.6%
Diffuse (Type II)	19.8%
Proliferative (Type III)	7.4%
Occlusive (Type IV)	1.2%

\* Composite of cardiac death, non-fatal Myocardial Infarction(MI), Clinically Driven Target Vessel Revascularization (CD-TVR)

\*\* Data shown as mean ± SD

## ISR distribution by stent type



## Acute and 6-month angiographic results\*

	Pre-procedure	Post-procedure	6-month
<b>In-stent</b>			
Reference vessel diameter	2.84 ± 0.39 mm	2.86 ± 0.39 mm	2.82 ± 0.38 mm
Minimum lumen diameter	0.91 ± 0.43 mm	2.18 ± 0.39 mm	2.08 ± 0.41 mm
<b>LLL</b>			<b>0.07 ± 0.31 mm</b>
Diameter stenosis	68.1 ± 13.8%	23.9 ± 9.8 %	25.9 ± 11.7 %
<b>In-segment</b>			
LLL			0.02 ± 0.32 mm

Primary endpoint

## 6-month clinical results

	All	BMS	DES
<b>MACE</b>	5 (6.5%)	2	3
Cardiac death	-	-	-
Non-fatal MI	1 (1.3%)	1	-
CD-TVR	4 (5.2%)	1	3
CD-TLR	3 (3.9%)	-	2

## 12-month clinical results

	All	BMS	DES
<b>MACE</b>	9 (11.8%)	2	7
Cardiac death	-	-	-
Non-fatal MI	1 (1.3%)	1	-
CD-TVR	8 (10.5%)	1	7
CD-TLR	7 (9.2%)	1	6

\*Data shown as mean ±SD

## Principal investigator

Prof. Ch. Hehrlein, University Medical Center Freiburg, Germany

Reference: Hehrlein C, et al. Twelve-month results of a paclitaxel releasing balloon in patients presenting with in-stent restenosis First-in-Man (PEPPER) trial. Cardiovasc Revasc Med. 2012 Sep; 13 (5): 260-264.