


Vascular Intervention // **Coronary**
Dual Lumen OTW Microcatheter

ReCross

A unique solution for guide wire redirection and support



Two separate OTW lumens and three exit ports allowing for exchange of both guide wires at any time¹



360° guide wire redirection options from two OTW lumens¹



Low profile, 6F trapping compatible¹



BIOTRONIK
excellence for life

ReCross

Dual Lumen OTW Microcatheter


Indicated to support guide wires during access of coronary vasculature.^a

Vascular
Intervention
Coronary



Uniquely versatile dual lumen catheter

The crossability of ReCross is comparable to best in class antegrade single lumen microcatheters. Two completely separate over the wire (OTW) lumens and a total of three exit ports, provide unique possibilities to resolve complex PCI.



Technical Data

Guide wire

Guide wire size	0.014"
Shaft design	Reinforced shaft
Lesion entry profile	1.5F
Exit markers	95 cm and 105 cm
Coating	Hydrophilic
Shaft profile	Oval shaft 2.3F x 3.3F
Guide wire compatibility	0.014"

Ordering Information

Ordering number	IMDS article number	Guiding catheter compatibility	Usable length
451671	RC1414025	5F	140 cm

^a Indication as per IFU.

1. IMDS data on file.

ReCross is a product of IMDS. Distributed by BIOTRONIK in selected countries.

RECCROSS

Manufacturer:
IMDS Operations B.V.
Ceintuurbaan Noord 150
9301 NZ Roden,
The Netherlands
Tel +31 (0)50 8200230
info@imds.nl

Distributor:
BIOTRONIK AG
Ackerstrasse 6
8180 Bülach, Switzerland
Tel +41 (0) 44 8645111
Fax +41 (0) 44 8645005
info.vi@biotronik.com
www.biotronik.com



© 2023 BIOTRONIK AG – All rights reserved.
Specifications are subject to modification, revision and improvement.

BIOTRONIK
excellence for life