



# UNIQUE SOLUTIONS FOR EVERYDAY COMPLEXITY

**BIOTRONIK Vascular Intervention** Product Catalog

For Professional Use Only. Product details presented here are intended for healthcare professionals in the United States.

# TABLE OF CONTENTS

## CORONARY VASCULAR INTERVENTION

<b>Orsiro Mission DES Sirolimus-Eluting Coronary Stent System</b> .....	<b>1</b>	<b>TrapIt Trapping Balloon Catheter</b> .....	<b>11</b>
Orsiro Mission 2.25 mm .....	1	TrapIt 90 cm .....	11
Orsiro Mission 2.50 mm .....	1	TrapIt 100 cm .....	11
Orsiro Mission 2.75 mm .....	2	<b>NHancer Rx Dual Lumen Microcatheter</b> .....	<b>12</b>
Orsiro Mission 3.00 mm .....	2	NHancer Rx.....	12
Orsiro Mission 3.50 mm .....	3	<b>ReCross Dual Lumen OTW Microcatheter</b> ....	<b>12</b>
Orsiro Mission 4.00 mm .....	4	ReCross .....	12
<b>PK Papyrus Covered Coronary Stent System</b> .....	<b>4</b>	<b>Micro Rx Rapid Exchange Microcatheter</b> .....	<b>12</b>
PK Papyrus 2.50 mm.....	4	Micro Rx.....	12
PK Papyrus 3.00 mm.....	4		
PK Papyrus 3.50 mm.....	5		
PK Papyrus 4.00 mm.....	5		
PK Papyrus 4.50 mm.....	5		
PK Papyrus 5.00 mm.....	6		
<b>Pantera Pro Workhorse Coronary Balloon</b> .....	<b>6</b>		
Pantera Pro 6 mm.....	6		
Pantera Pro 10 mm.....	6		
Pantera Pro 15 mm.....	7		
Pantera Pro 20 mm.....	7		
Pantera Pro 25 mm.....	7		
Pantera Pro 30 mm.....	8		
<b>Pantera LEO Non-Compliant Coronary Balloon</b> .....	<b>8</b>		
Pantera LEO 8 mm.....	8		
Pantera LEO 12 mm.....	9		
Pantera LEO 15 mm.....	9		
Pantera LEO 20 mm.....	10		
Pantera LEO 30 mm.....	11		








# PERIPHERAL VASCULAR INTERVENTION

<b>Oscar Peripheral Multifunctional Catheter System</b> .....	<b>13</b>	<b>Passeo-35 Xeo Peripheral PTA Balloon Catheter</b> .....	<b>32</b>
Oscar 0.014"/4F Multifunctional Catheter System .....	13	Passeo-35 Xeo 90 cm Catheter Length .....	32
Single Oscar PTA Balloons 0.014"/4F .....	13	Passeo-35 Xeo 130 cm Catheter Length .....	35
Oscar 0.018"/6F Multifunctional Catheter System .....	13	Passeo-35 Xeo 170 cm Catheter Length .....	39
Single Oscar PTA Balloons 0.018"/6F.....	14		
<b>Pulsar-18 T3 Peripheral Self-Expanding Stent System</b> .....	<b>14</b>	<b>Passeo-35 HP Peripheral High Pressure PTA Balloon Catheter</b> .....	<b>41</b>
Pulsar-18 T3 90 cm Catheter Length .....	14	Passeo-35 HP 40 cm Catheter Length .....	41
Pulsar-18 T3 135 cm Catheter Length .....	17	Passeo-35 HP 75 cm Catheter Length .....	42
<b>Astron Peripheral Self-Expanding Stent System</b> .....	<b>20</b>	<b>Fortress Peripheral Reinforced Introducer Sheath System</b> .....	<b>43</b>
Astron 72 cm Catheter Length .....	20	Fortress Introducer Sheath .....	43
Astron 130 cm Catheter Length .....	21		
<b>Passeo-14 Peripheral PTA Balloon Catheter</b> .....	<b>21</b>		
Passeo-14 90 cm Catheter Length .....	21		
Passeo-14 120 cm Catheter Length .....	23		
Passeo-14 150 cm Catheter Length .....	23		
<b>Passeo-18 Peripheral PTA Balloon Catheter</b> .....	<b>25</b>		
Passeo-18 90 cm Catheter Length .....	25		
Passeo-18 130 cm Catheter Length .....	29		
Passeo-18 150 cm Catheter Length .....	32		








# CORONARY VASCULAR INTERVENTION

## Orsiro Mission DES Sirolimus-Eluting Coronary Stent System



### Orsiro Mission 2.25 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453925	2.25	9	10	16	C1874	 010764013045581817453925
453931	2.25	13	10	16	C1874	 010764013045587017453931
453937	2.25	15	10	16	C1874	 010764013045593117453937
453943	2.25	18	10	16	C1874	 010764013045599317453943
453949	2.25	22	10	16	C1874	 010764013045605117453949
453955	2.25	26	10	16	C1874	 010764013045611217453955
453961	2.25	30	10	16	C1874	 010764013045617417453961










### Orsiro Mission 2.50 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453926	2.5	9	10	16	C1874	 010764013045582517453926
453932	2.5	13	10	16	C1874	 010764013045588717453932
453938	2.5	15	10	16	C1874	 010764013045594817453938
453944	2.5	18	10	16	C1874	 010764013045600617453944
453950	2.5	22	10	16	C1874	 010764013045606817453950
453956	2.5	26	10	16	C1874	 010764013045612917453956
453962	2.5	30	10	16	C1874	 010764013045618117453962




### Orsiro Mission 2.50 mm (Continued)

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453968	2.5	35	10	16	C1874	 010764013045624217453968
453974	2.5	40	10	16	C1874	 010764013045630317453974







### Orsiro Mission 2.75 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453927	2.75	9	10	16	C1874	 010764013045583217453927
453933	2.75	13	10	16	C1874	 010764013045589417453933
453939	2.75	15	10	16	C1874	 010764013045595517453939
453945	2.75	18	10	16	C1874	 010764013045601317453945
453951	2.75	22	10	16	C1874	 010764013045607517453951
453957	2.75	26	10	16	C1874	 010764013045613617453957
453963	2.75	30	10	16	C1874	 010764013045619817453963
453969	2.75	35	10	16	C1874	 010764013045625917453969
453975	2.75	40	10	16	C1874	 010764013045631017453975










### Orsiro Mission 3.00 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453928	3	9	10	16	C1874	 010764013045584917453928
453934	3	13	10	16	C1874	 010764013045590017453934
453940	3	15	10	16	C1874	 010764013045596217453940










### Orsiro Mission 3.00 mm (Continued)

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453946	3	18	10	16	C1874	 010764013045602017453946
453952	3	22	10	16	C1874	 010764013045608217453952
453958	3	26	10	16	C1874	 010764013045614317453958
453964	3	30	10	16	C1874	 010764013045620417453964
453970	3	35	10	16	C1874	 010764013045626617453970
453976	3	40	10	16	C1874	 010764013045632717453976

### Orsiro Mission 3.50 mm



ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453929	3.5	9	10	16	C1874	 010764013045585617453929
453935	3.5	13	10	16	C1874	 010764013045591717453935
453941	3.5	15	10	16	C1874	 010764013045597917453941
453947	3.5	18	10	16	C1874	 010764013045603717453947
453953	3.5	22	10	16	C1874	 010764013045609917453953
453959	3.5	26	10	16	C1874	 010764013045615017453959
453965	3.5	30	10	16	C1874	 010764013045621117453965
453971	3.5	35	10	16	C1874	 010764013045627317453971
453977	3.5	40	10	16	C1874	 010764013045633417453977

## Orsiro Mission 4.00 mm


ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
453930	4	9	10	16	C1874	 010764013045586317453930
453936	4	13	10	16	C1874	 010764013045592417453936
453942	4	15	10	16	C1874	 010764013045598617453942
453948	4	18	10	16	C1874	 010764013045604417453948
453954	4	22	10	16	C1874	 010764013045610517453954
453960	4	26	10	16	C1874	 010764013045616717453960
453966	4	30	10	16	C1874	 010764013045622817453966
453972	4	35	10	16	C1874	 010764013045628017453972
453978	4	40	10	16	C1874	 010764013045634117453978

## PK Papyrus Covered Coronary Stent System



### PK Papyrus 2.50 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434887	2.50	15	8	16	C1874	 010764013044715817434887
434893	2.50	20	8	16	C1874	 010764013044721917434893




### PK Papyrus 3.00 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434888	3.00	15	8	16	C1874	 010764013044716517434888




### PK Papyrus 3.00 mm (Continued)

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434894	3.00	20	8	16	C1874	 010764013044722617434894
434899	3.00	26	8	16	C1874	 01076401304472117434899




### PK Papyrus 3.50 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434889	3.50	15	8	16	C1874	 010764013044717217434889
434895	3.50	20	8	16	C1874	 010764013044723317434895
434900	3.50	26	8	16	C1874	 010764013044728817434900

### PK Papyrus 4.00 mm




ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434890	4.00	15	7	16	C1874	 010764013044718917434890
434896	4.00	20	7	16	C1874	 010764013044724017434896
434901	4.00	26	7	16	C1874	 010764013044729517434901

### PK Papyrus 4.50 mm

ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434891	4.50	15	7	14	C1874	 010764013044719617434891
434897	4.50	20	7	14	C1874	 010764013044725717434897
434902	4.50	26	7	14	C1874	 010764013044730117434902





## PK Papyrus 5.00 mm








ORDER NUMBER	STENT DIAMETER (MM)	STENT LENGTH (MM)	NOMINAL PRESSURE (ATM)	RATED BURST PRESSURE (ATM)	C-CODE	BARCODE (GTIN)
434892	5.00	15	7	14	C1874	 010764013044720217434892
434898	5.00	20	7	14	C1874	 010764013044726417434898
434903	5.00	26	7	14	C1874	 010764013044731817434903

## Pantera Pro Workhorse Coronary Balloon








### Pantera Pro 6 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393289	1.25	6	140	C1725	 010764013043693017393289
393290	1.50	6	140	C1725	 010764013043694717393290








### Pantera Pro 10 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393291	1.25	10	140	C1725	 010764013043695417393291
393292	1.50	10	140	C1725	 010764013043696117393292
393293	2.00	10	140	C1725	 010764013043697817393293
393294	2.50	10	140	C1725	 010764013043698517393294
393295	3.00	10	140	C1725	 010764013043699217393295
393296	3.50	10	140	C1725	 010764013043700517393296
393297	4.00	10	140	C1725	 010764013043701217393297




## Pantera Pro 15 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393298	1.25	15	140	C1725	 010764013043702917393298
393299	1.50	15	140	C1725	 010764013043703617393299
393300	2.00	15	140	C1725	 010764013043704317393300
393301	2.50	15	140	C1725	 010764013043705017393301
393302	3.00	15	140	C1725	 010764013043706717393302
393303	3.50	15	140	C1725	 010764013043707417393303
393304	4.00	15	140	C1725	 010764013043708117393304



## Pantera Pro 20 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393305	1.25	20	140	C1725	 010764013043709817393305
393306	1.50	20	140	C1725	 010764013043710417393306
393307	2.00	20	140	C1725	 010764013043711117393307
393308	2.50	20	140	C1725	 010764013043712817393308
393309	3.00	20	140	C1725	 010764013043713517393309
393310	3.50	20	140	C1725	 010764013043714217393310
393311	4.00	20	140	C1725	 010764013043715917393311






## Pantera Pro 25 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393312	2.00	25	140	C1725	 010764013043716617393312
393313	2.50	25	140	C1725	 010764013043717317393313
393314	3.00	25	140	C1725	 010764013043718017393314

## Pantera Pro 25mm (Continued)









ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393315	3.50	25	140	C1725	 010764013043719717393315
393316	4.00	25	140	C1725	 010764013043720317393316

## Pantera Pro 30 mm




ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
393317	2.00	30	140	C1725	 010764013043721017393317
393318	2.50	30	140	C1725	 010764013043722717393318
393319	3.00	30	140	C1725	 010764013043723417393319
393320	3.50	30	140	C1725	 010764013043724117393320
393321	4.00	30	140	C1725	 010764013043725817393321

## Pantera LEO Non-Compliant Coronary Balloon

### Pantera LEO 8 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
366991	2.00	8	145	C1725	 010764013041549217366991
366992	2.25	8	145	C1725	 010764013041550817366992
366993	2.50	8	145	C1725	 010764013041551517366993
366994	2.75	8	145	C1725	 010764013041552217366994
366995	3.00	8	145	C1725	 010764013041553917366995
366996	3.25	8	145	C1725	 010764013041554617366996
366997	3.50	8	145	C1725	 010764013041555317366997
366998	3.75	8	145	C1725	 010764013041556017366998




## Pantera LEO 8 mm (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
366999	4.00	8	145	C1725	 010764013041557717366999
367000	4.50	8	145	C1725	 010764013041558417367000
367001	5.00	8	145	C1725	 010764013041559117367001









## Pantera LEO 12 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
367002	2.00	12	145	C1725	 010764013041560717367002
367003	2.25	12	145	C1725	 010764013041561417367003
367004	2.50	12	145	C1725	 010764013041562117367004
367005	2.75	12	145	C1725	 010764013041563817367005
367006	3.00	12	145	C1725	 010764013041564517367006
367007	3.25	12	145	C1725	 010764013041565217367007
367008	3.50	12	145	C1725	 010764013041566917367008
367009	3.75	12	145	C1725	 010764013041567617367009
367010	4.00	12	145	C1725	 010764013041568317367010
367011	4.50	12	145	C1725	 010764013041569017367011
367012	5.00	12	145	C1725	 010764013041570617367012












## Pantera LEO 15 mm

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
367013	2.00	15	145	C1725	 010764013041571317367013
367014	2.25	15	145	C1725	 010764013041572017367014
367015	2.50	15	145	C1725	 010764013041573717367015










## Pantera LEO 15 mm (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
367016	2.75	15	145	C1725	 010764013041574417367016
367017	3.00	15	145	C1725	 010764013041575117367017
367018	3.25	15	145	C1725	 010764013041576817367018
367019	3.50	15	145	C1725	 010764013041577517367019
367020	3.75	15	145	C1725	 010764013041578217367020
367021	4.00	15	145	C1725	 010764013041579917367021
367022	4.50	15	145	C1725	 010764013041580517367022
367023	5.00	15	145	C1725	 010764013041581217367023

## Pantera LEO 20 mm


ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
367024	2.00	20	145	C1725	 010764013041582917367024
367025	2.25	20	145	C1725	 010764013041583617367025
367026	2.50	20	145	C1725	 010764013041584317367026
367027	2.75	20	145	C1725	 010764013041585017367027
367028	3.00	20	145	C1725	 010764013041586717367028
367029	3.25	20	145	C1725	 010764013041587417367029
367030	3.50	20	145	C1725	 010764013041588117367030
367031	3.75	20	145	C1725	 010764013041589817367031
367032	4.00	20	145	C1725	 010764013041590417367032
367033	4.50	20	145	C1725	 010764013041591117367033
367034	5.00	20	145	C1725	 010764013041592817367034

## Pantera LEO 30 mm


ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	EFFECTIVE LENGTH (CM)	C-CODE	BARCODE (GTIN)
367035	2.00	30	145	C1725	 010764013041593517367035
367036	2.25	30	145	C1725	 010764013041594217367036
367037	2.50	30	145	C1725	 010764013041595917367037
367038	2.75	30	145	C1725	 010764013041596617367038
367039	3.00	30	145	C1725	 010764013041597317367039
367040	3.25	30	145	C1725	 010764013041598017367040
367041	3.50	30	145	C1725	 010764013041599717367041
367042	3.75	30	145	C1725	 010764013041600017367042
367043	4.00	30	145	C1725	 010764013041601717367043
367044	4.50	30	145	C1725	 010764013041602417367044
367045	5.00	30	145	C1725	 010764013041603117367045

## TrapIt Trapping Balloon Catheter

### TrapIt 90 cm


ORDER NUMBER	DESCRIPTION	C-CODE	BARCODE (GTIN)
451667	Trapping balloon	C1725	 010871848174026217451667

### TrapIt 100 cm

ORDER NUMBER	DESCRIPTION	C-CODE	BARCODE (GTIN)
451668	Trapping balloon	C1725	 010871848174027917451668


## NHancer Rx Dual Lumen Microcatheter

### NHancer Rx

ORDER NUMBER	LESION ENTRY PROFILE	DISTAL PROFILE	USABLE LENGTH (CM)	GUIDE WIRE COMPATIBILITY (IN)	C-CODE	BARCODE (GTIN)
451669	1.5F	Oval low profile design 2.3F x 3.3F	135	0.014"	C1725	 010871848174023117451669


## ReCross Dual Lumen OTW Microcatheter

### ReCross

ORDER NUMBER	LESION ENTRY PROFILE	SHAFT PROFILE	USABLE LENGTH (CM)	GUIDE WIRE COMPATIBILITY (IN)	C-CODE	BARCODE (GTIN)
451671	1.5F	Oval shaft 2.3F x 3.3F	140	0.014"	C1725	 010871848174030917451671

## Micro Rx Rapid Exchange Microcatheter





### Micro Rx

ORDER NUMBER	LESION ENTRY PROFILE	USABLE LENGTH (CM)	GUIDE WIRE COMPATIBILITY (IN)	C-CODE	BARCODE (GTIN)
483054	1.3F	135	0.014"	C1725	 010871848174035417483054

# PERIPHERAL VASCULAR INTERVENTION





## Oscar Peripheral Multifunctional Catheter System

### Oscar 0.014"/4F Multifunctional Catheter System




ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
454332	2.0	20-180	16	127	4	0.014"	C7531	 010764013045559717454332
454333	2.5	20-180	16	127	4	0.014"	C7531	 010764013045580117454333
454334	3.0	20-180	16	127	4	0.014"	C7531	 010764013045579517454334
454335	3.5	20-180	14	127	4	0.014"	C7531	 010764013045578817454335

### Single Oscar PTA Balloons 0.014"/4F

\*Only to be used with the matching Oscar Multifunctional Support Catheter System




ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
454343	2.0	20-180	16	127	4	0.014"	C7531	 010764013045570217454343
454344	2.5	20-180	16	127	4	0.014"	C7531	 010764013045569617454344
454345	3.0	20-180	16	127	4	0.014"	C7531	 010764013045568917454345
454346	3.5	20-180	14	127	4	0.014"	C7531	 010764013045567217454346

### Oscar 0.018"/6F Multifunctional Catheter System

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
454336	4.0	20-180	14	128	6	0.018"	C7531	 010764013045577117454336
454337	5.0	20-180	14	128	6	0.018"	C7531	 010764013045576417454337
454338	6.0	20-180	14	128	6	0.018"	C7531	 010764013045575717454338




## Oscar 0.018"/6F Multifunctional Catheter System (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
454339	4.0	20-180	14	80	6	0.018"	C7531	 010764013045574017454339
454340	5.0	20-180	14	80	6	0.018"	C7531	 010764013045573317454340
454341	6.0	20-180	14	80	6	0.018"	C7531	 010764013045572617454341



## Single Oscar PTA Balloons 0.018"/6F

\*Only to be used with the matching Oscar Multifunctional Support Catheter System

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
454347	4.0	20-180	14	128	6	0.018"	C7531	 010764013045566517454347
454348	5.0	20-180	14	128	6	0.018"	C7531	 010764013045565817454348
454349	6.0	20-180	14	128	6	0.018"	C7531	 010764013045564117454349
454350	4.0	20-180	14	80	6	0.018"	C7531	 010764013045563417454350
454351	5.0	20-180	14	80	6	0.018"	C7531	 010764013045562717454351
454352	6.0	20-180	14	80	6	0.018"	C7531	 010764013045561017454352

## Pulsar-18 T3 Peripheral Self-Expanding Stent System


### Pulsar-18 T3 90 cm Catheter Length

ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430437	4	20	90	4	0.018"	C1876	 010764013044635917430437
430438	4	30	90	4	0.018"	C1876	 010764013044636617430438





Pulsar-18 T3 90 cm Catheter Length (Continued)

ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430439	4	40	90	4	0.018"	C1876	 010764013044637317430439
430440	4	60	90	4	0.018"	C1876	 010764013044638017430440
430441	4	80	90	4	0.018"	C1876	 010764013044639717430441
430442	4	100	90	4	0.018"	C1876	 010764013044640317430442
430443	4	120	90	4	0.018"	C1876	 010764013044641017430443
430444	4	150	90	4	0.018"	C1876	 010764013044642717430444
430445	4	170	90	4	0.018"	C1876	 010764013044643417430445
430446	4	200	90	4	0.018"	C1876	 010764013044644117430446
430447	5	20	90	4	0.018"	C1876	 010764013044645817430447
430448	5	30	90	4	0.018"	C1876	 010764013044646517430448
430449	5	40	90	4	0.018"	C1876	 010764013044647217430449
430450	5	60	90	4	0.018"	C1876	 010764013044648917430450
430451	5	80	90	4	0.018"	C1876	 010764013044649617430451
430452	5	100	90	4	0.018"	C1876	 010764013044650217430452
430453	5	120	90	4	0.018"	C1876	 010764013044651917430453
430454	5	150	90	4	0.018"	C1876	 010764013044652617430454
430455	5	170	90	4	0.018"	C1876	 010764013044653317430455











Pulsar-18 T3 90 cm Catheter Length (Continued)

ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430456	5	200	90	4	0.018"	C1876	 010764013044654017430456
430457	6	20	90	4	0.018"	C1876	 010764013044655717430457
430458	6	30	90	4	0.018"	C1876	 010764013044656417430458
430459	6	40	90	4	0.018"	C1876	 010764013044657117430459
430460	6	60	90	4	0.018"	C1876	 010764013044658817430460
430461	6	80	90	4	0.018"	C1876	 010764013044659517430461
430462	6	100	90	4	0.018"	C1876	 010764013044660117430462
430463	6	120	90	4	0.018"	C1876	 010764013044661817430463
430464	6	150	90	4	0.018"	C1876	 010764013044662517430464
430465	6	170	90	4	0.018"	C1876	 010764013044663217430465
430466	6	200	90	4	0.018"	C1876	 010764013044664917430466
430467	7	20	90	4	0.018"	C1876	 010764013044665617430467
430468	7	30	90	4	0.018"	C1876	 010764013044666317430468
430469	7	40	90	4	0.018"	C1876	 010764013044667017430469
430470	7	60	90	4	0.018"	C1876	 010764013044668717430470
430471	7	80	90	4	0.018"	C1876	 010764013044669417430471
430472	7	100	90	4	0.018"	C1876	 010764013044670017430472

**Pulsar-18 T3 90 cm Catheter Length (Continued)**

ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430473	7	120	90	4	0.018"	C1876	 010764013044671717430473
430474	7	150	90	4	0.018"	C1876	 010764013044672417430474
430475	7	170	90	4	0.018"	C1876	 010764013044673117430475
430476	7	200	90	4	0.018"	C1876	 010764013044674817430476














**Pulsar-18 T3 135 cm Catheter Length**

ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430477	4	20	135	4	0.018"	C1876	 010764013044675517430477
430478	4	30	135	4	0.018"	C1876	 010764013044676217430478
430479	4	40	135	4	0.018"	C1876	 010764013044677917430479
430480	4	60	135	4	0.018"	C1876	 010764013044678617430480
430481	4	80	135	4	0.018"	C1876	 010764013044679317430481
430482	4	100	135	4	0.018"	C1876	 010764013044680917430482
430483	4	120	135	4	0.018"	C1876	 010764013044681617430483
430484	4	150	135	4	0.018"	C1876	 010764013044682317430484
430485	4	170	135	4	0.018"	C1876	 010764013044683017430485
430486	4	200	135	4	0.018"	C1876	 010764013044684717430486

**Pulsar-18 T3 135 cm Catheter Length (Continued)**
















ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430487	5	20	135	4	0.018"	C1876	 010764013044685417430487
430488	5	30	135	4	0.018"	C1876	 010764013044686117430488
430489	5	40	135	4	0.018"	C1876	 010764013044687817430489
430490	5	60	135	4	0.018"	C1876	 010764013044688517430490
430491	5	80	135	4	0.018"	C1876	 010764013044689217430491
430492	5	100	135	4	0.018"	C1876	 010764013044690817430492
430493	5	120	135	4	0.018"	C1876	 010764013044691517430493
430494	5	150	135	4	0.018"	C1876	 010764013044692217430494
430495	5	170	135	4	0.018"	C1876	 010764013044693917430495
430496	5	200	135	4	0.018"	C1876	 010764013044694617430496
430497	6	20	135	4	0.018"	C1876	 010764013044695317430497
430498	6	30	135	4	0.018"	C1876	 010764013044696017430498
430499	6	40	135	4	0.018"	C1876	 010764013044697717430499
430500	6	60	135	4	0.018"	C1876	 010764013044698417430500
430501	6	80	135	4	0.018"	C1876	 010764013044699117430501
430502	6	100	135	4	0.018"	C1876	 010764013044700417430502
430503	6	120	135	4	0.018"	C1876	 010764013044701117430503

**Pulsar-18 T3 135 cm Catheter Length (Continued)**













ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
430504	6	150	135	4	0.018"	C1876	 010764013044702817430504
430505	6	170	135	4	0.018"	C1876	 010764013044703517430505
430506	6	200	135	4	0.018"	C1876	 010764013044704217430506
430507	7	20	135	4	0.018"	C1876	 010764013044705917430507
430508	7	30	135	4	0.018"	C1876	 010764013044706617430508
430509	7	40	135	4	0.018"	C1876	 010764013044707317430509
430510	7	60	135	4	0.018"	C1876	 010764013044708017430510
430511	7	80	135	4	0.018"	C1876	 010764013044709717430511
430512	7	100	135	4	0.018"	C1876	 010764013044710317430512
430513	7	120	135	4	0.018"	C1876	 010764013044711017430513
430514	7	150	135	4	0.018"	C1876	 010764013044712717430514
430515	7	170	135	4	0.018"	C1876	 010764013044713417430515
430516	7	200	135	4	0.018"	C1876	 010764013044714117430516

# Astron Peripheral Self-Expanding Stent System

## Astron 72 cm Catheter Length



ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
364645	7	30	72	6	0.035"	C1876	 010764013041220017364645
364646	7	40	72	6	0.035"	C1876	 010764013041221717364646
364647	7	60	72	6	0.035"	C1876	 010764013041222417364647
364648	7	80	72	6	0.035"	C1876	 010764013041223117364648
364649	8	30	72	6	0.035"	C1876	 010764013041224817364649
364650	8	40	72	6	0.035"	C1876	 010764013041225517364650
364651	8	60	72	6	0.035"	C1876	 010764013041226217364651
364652	8	80	72	6	0.035"	C1876	 010764013041227917364652
364653	9	30	72	6	0.035"	C1876	 010764013041228617364653
364654	9	40	72	6	0.035"	C1876	 010764013041229317364654
364655	9	60	72	6	0.035"	C1876	 010764013041230917364655
364656	9	80	72	6	0.035"	C1876	 010764013041231617364656
364657	10	40	72	6	0.035"	C1876	 010764013041232317364657
364658	10	60	72	6	0.035"	C1876	 010764013041233017364658
364659	10	80	72	6	0.035"	C1876	 010764013041234717364659

## Astron 130 cm Catheter Length

ORDER NUMBER	STENT DIAMETER (MM)	NOMINAL STENT LENGTH (MM)	USABLE LENGTH (CM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
364660	7	30	130	6	0.035"	C1876	 010764013041235417364660
364661	7	40	130	6	0.035"	C1876	 010764013041236117364661
364662	7	60	130	6	0.035"	C1876	 010764013041237817364662
364663	7	80	130	6	0.035"	C1876	 010764013041238517364663
364664	8	30	130	6	0.035"	C1876	 010764013041239217364664
364665	8	40	130	6	0.035"	C1876	 010764013041240817364665
364666	8	60	130	6	0.035"	C1876	 010764013041241517364666
364667	8	80	130	6	0.035"	C1876	 010764013041242217364667
364668	9	30	130	6	0.035"	C1876	 010764013041243917364668
364669	9	40	130	6	0.035"	C1876	 010764013041244617364669
364670	9	60	130	6	0.035"	C1876	 010764013041245317364670
364671	9	80	130	6	0.035"	C1876	 010764013041246017364671














## Passeo-14 Peripheral PTA Balloon Catheter

### Passeo-14 90 cm Catheter Length










ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
380279	2.5	40	14	4	0.014"	C1725	 010764013043315117380279
380280	3.0	40	14	4	0.014"	C1725	 010764013043316817380280









### Passeo-14 90 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
380281	3.5	40	14	4	0.014"	C1725	 010764013043317517380281
380282	4.0	40	14	4	0.014"	C1725	 010764013043318217380282
380285	2.5	70	14	4	0.014"	C1725	 010764013043321217380285
380286	3.0	70	14	4	0.014"	C1725	 010764013043322917380286
380287	3.5	70	14	4	0.014"	C1725	 010764013043323617380287
380288	4.0	70	14	4	0.014"	C1725	 010764013043324317380288
380291	2.5	100	14	4	0.014"	C1725	 010764013043327417380291
380292	3.0	100	14	4	0.014"	C1725	 010764013043328117380292
380293	3.5	100	14	4	0.014"	C1725	 010764013043329817380293
380294	4.0	100	14	4	0.014"	C1725	 010764013043330417380294
380297	2.5	140	14	4	0.014"	C1725	 010764013043333517380297
380298	3.0	140	14	4	0.014"	C1725	 010764013043334217380298
380299	3.5	140	14	4	0.014"	C1725	 010764013043335917380299
380300	4.0	140	14	4	0.014"	C1725	 010764013043336617380300
380303	2.5	180	14	4	0.014"	C1725	 010764013043339717380303
380304	3.0	180	14	4	0.014"	C1725	 010764013043340317380304
380309	2.5	220	14	4	0.014"	C1725	 010764013043345817380309
380310	3.0	220	14	4	0.014"	C1725	 010764013043346517380310




















### Passeo-14 120 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
380271	1.5	20	14	4	0.014"	C1725	 010764013043307617380271
380277	1.5	40	14	4	0.014"	C1725	 010764013043313717380277
380278	2.0	40	14	4	0.014"	C1725	 010764013043314417380278
380283	1.5	70	14	4	0.014"	C1725	 010764013043319917380283
380284	2.0	70	14	4	0.014"	C1725	 010764013043320517380284
380290	2.0	100	14	4	0.014"	C1725	 010764013043326717380290
380296	2.0	140	14	4	0.014"	C1725	 010764013043332817380296
380302	2.0	180	14	4	0.014"	C1725	 010764013043338017380302
380308	2.0	220	14	4	0.014"	C1725	 010764013043344117380308





### Passeo-14 150 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
380313	1.5	20	14	4	0.014"	C1725	 010764013043349617380313
380319	1.5	40	14	4	0.014"	C1725	 010764013043355717380319
380320	2.0	40	14	4	0.014"	C1725	 010764013043356417380320
380321	2.5	40	14	4	0.014"	C1725	 010764013043357117380321
380322	3.0	40	14	4	0.014"	C1725	 010764013043358817380322
380323	3.5	40	14	4	0.014"	C1725	 010764013043359517380323

## Passeo-14 150 cm Catheter Length (Continued)












ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
380324	4.0	40	14	4	0.014"	C1725	 010764013043360117380324
380325	1.5	70	14	4	0.014"	C1725	 010764013043361817380325
380326	2.0	70	14	4	0.014"	C1725	 010764013043362517380326
380327	2.5	70	14	4	0.014"	C1725	 010764013043363217380327
380328	3.0	70	14	4	0.014"	C1725	 010764013043364917380328
380329	3.5	70	14	4	0.014"	C1725	 010764013043365617380329
380330	4.0	70	14	4	0.014"	C1725	 010764013043366317380330
380332	2.0	100	14	4	0.014"	C1725	 010764013043368717380332
380333	2.5	100	14	4	0.014"	C1725	 010764013043369417380333
380334	3.0	100	14	4	0.014"	C1725	 010764013043370017380334
380335	3.5	100	14	4	0.014"	C1725	 010764013043371717380335
380336	4.0	100	14	4	0.014"	C1725	 010764013043372417380336
380338	2.0	140	14	4	0.014"	C1725	 010764013043374817380338
380339	2.5	140	14	4	0.014"	C1725	 010764013043375517380339
380340	3.0	140	14	4	0.014"	C1725	 010764013043376217380340
380341	3.5	140	14	4	0.014"	C1725	 010764013043377917380341
380342	4.0	140	14	4	0.014"	C1725	 010764013043378617380342
380344	2.0	180	14	4	0.014"	C1725	 010764013043380917380344
380345	2.5	180	14	4	0.014"	C1725	 010764013043381617380345

## Passeo-14 150 cm Catheter Length (Continued)


















ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
380346	3.0	180	14	4	0.014"	C1725	 010764013043382317380346
380350	2.0	220	14	4	0.014"	C1725	 010764013043386117380350
380351	2.5	220	14	4	0.014"	C1725	 010764013043387817380351
380352	3.0	220	14	4	0.014"	C1725	 010764013043388517380352

## Passeo-18 Peripheral PTA Balloon Catheter




















### Passeo-18 90 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
366098	2.0	20	15	4	0.018"	C1725	 010764013041440217366098
357451	2.5	20	15	4	0.018"	C1725	 010764011955862217357451
357452	3.0	20	15	4	0.018"	C1725	 010764011955863917357452
357453	3.5	20	15	4	0.018"	C1725	 010764011955864617357453
357454	4.0	20	15	4	0.018"	C1725	 010764011955865317357454
357455	5.0	20	15	4	0.018"	C1725	 010764011955866017357455
357456	6.0	20	12	4	0.018"	C1725	 010764011955867717357456
357457	7.0	20	12	4	0.018"	C1725	 010764011955868417357457
366099	2.0	40	15	4	0.018"	C1725	 010764013041441917366099
357458	2.5	40	15	4	0.018"	C1725	 010764011955869117357458
357459	3.0	40	15	4	0.018"	C1725	 010764011955870717357459
















Passeo-18 90 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
357460	3.5	40	15	4	0.018"	C1725	 010764011955871417357460
357461	4.0	40	15	4	0.018"	C1725	 010764011955872117357461
357462	5.0	40	15	4	0.018"	C1725	 010764011955873817357462
357463	6.0	40	12	4	0.018"	C1725	 010764011955874517357463
357464	7.0	40	12	4	0.018"	C1725	 010764011955875217357464
366100	2.0	60	15	4	0.018"	C1725	 010764013041442617366100
366101	2.5	60	15	4	0.018"	C1725	 010764013041443317366101
366102	3.0	60	15	4	0.018"	C1725	 010764013041444017366102
366103	3.5	60	15	4	0.018"	C1725	 010764013041445717366103
357465	4.0	60	15	4	0.018"	C1725	 010764011955876917357465
357466	5.0	60	15	4	0.018"	C1725	 010764011955877617357466
357467	6.0	60	12	4	0.018"	C1725	 010764011955878317357467
357468	7.0	60	12	4	0.018"	C1725	 010764011955879017357468
366104	2.0	80	15	4	0.018"	C1725	 010764013041446417366104
357469	2.5	80	15	4	0.018"	C1725	 010764011955880617357469
357470	3.0	80	15	4	0.018"	C1725	 010764011955881317357470
357471	3.5	80	15	4	0.018"	C1725	 010764011955882017357471
357472	4.0	80	15	4	0.018"	C1725	 010764011955883717357472
357473	5.0	80	15	4	0.018"	C1725	 010764011955884417357473




Passeo-18 90 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
357474	6.0	80	12	4	0.018"	C1725	 010764011955885117357474
357475	7.0	80	12	5	0.018"	C1725	 010764011955886817357475
366105	2.0	120	15	4	0.018"	C1725	 010764013041447117366105
357476	2.5	120	15	4	0.018"	C1725	 010764011955887517357476
357477	3.0	120	15	4	0.018"	C1725	 010764011955888217357477
357478	3.5	120	15	4	0.018"	C1725	 010764011955889917357478
357479	4.0	120	15	4	0.018"	C1725	 010764011955890517357479
357480	5.0	120	15	4	0.018"	C1725	 010764011955891217357480
357481	6.0	120	12	5	0.018"	C1725	 010764011955892917357481
357482	7.0	120	12	5	0.018"	C1725	 010764011955893617357482
366106	2.0	150	15	4	0.018"	C1725	 010764013041448817366106
366107	2.5	150	15	4	0.018"	C1725	 010764013041449517366107
366108	3.0	150	15	4	0.018"	C1725	 010764013041450117366108
366109	3.5	150	15	4	0.018"	C1725	 010764013041451817366109
366110	4.0	150	15	4	0.018"	C1725	 010764013041452517366110
366111	5.0	150	12	4	0.018"	C1725	 010764013041453217366111
366112	6.0	150	12	5	0.018"	C1725	 010764013041454917366112
366113	7.0	150	12	5	0.018"	C1725	 010764013041455617366113
366114	2.0	170	15	4	0.018"	C1725	 010764013041456317366114

Passeo-18 90 cm Catheter Length (Continued)




















ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
357483	2.5	170	15	4	0.018"	C1725	 010764011955894317357483
357484	3.0	170	15	4	0.018"	C1725	 010764011955895017357484
357485	3.5	170	15	4	0.018"	C1725	 010764011955896717357485
376272	4.0	170	13	4	0.018"	C1725	 010764013042979617376272
376273	5.0	170	13	4	0.018"	C1725	 010764013042980217376273
376274	6.0	170	12	5	0.018"	C1725	 010764013042981917376274
376275	7.0	170	12	5	0.018"	C1725	 010764013042982617376275
376276	2.0	200	14	4	0.018"	C1725	 010764013042983317376276
376277	2.5	200	14	4	0.018"	C1725	 010764013042984017376277
376278	3.0	200	14	4	0.018"	C1725	 010764013042985717376278
376279	3.5	200	14	4	0.018"	C1725	 010764013042986417376279
376280	4.0	200	13	4	0.018"	C1725	 010764013042987117376280
376281	5.0	200	13	4	0.018"	C1725	 010764013042988817376281
376282	6.0	200	12	5	0.018"	C1725	 010764013042989517376282
376283	7.0	200	12	5	0.018"	C1725	 010764013042990117376283

## Passeo-18 130 cm Catheter Length







ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
357486	2.5	20	15	4	0.018"	C1725	 010764011955897417357486
357487	3.0	20	15	4	0.018"	C1725	 010764011955898117357487
357488	3.5	20	15	4	0.018"	C1725	 010764011955899817357488
357489	4.0	20	15	4	0.018"	C1725	 010764011955900117357489
357490	5.0	20	15	4	0.018"	C1725	 010764011955901817357490
366116	6.0	20	12	4	0.018"	C1725	 010764013041458717366116
366117	7.0	20	12	4	0.018"	C1725	 010764013041459417366117
357491	2.5	40	15	4	0.018"	C1725	 010764011955902517357491
357492	3.0	40	15	4	0.018"	C1725	 010764011955903217357492
357493	3.5	40	15	4	0.018"	C1725	 010764011955904917357493
357494	4.0	40	15	4	0.018"	C1725	 010764011955905617357494
357495	5.0	40	15	4	0.018"	C1725	 010764011955906317357495
357496	6.0	40	12	4	0.018"	C1725	 010764011955907017357496
357497	7.0	40	12	4	0.018"	C1725	 010764011955908717357497
366120	2.5	60	15	4	0.018"	C1725	 010764013041462417366120
366121	3.0	60	15	4	0.018"	C1725	 010764013041463117366121
366122	3.5	60	15	4	0.018"	C1725	 010764013041464817366122
357498	4.0	60	15	4	0.018"	C1725	 010764011955909417357498
357499	5.0	60	15	4	0.018"	C1725	 010764011955910017357499







Passeo-18 130 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
357500	6.0	60	12	4	0.018"	C1725	 010764011955911717357500
357501	7.0	60	12	4	0.018"	C1725	 010764011955912417357501
357502	2.5	80	15	4	0.018"	C1725	 010764011955913117357502
357503	3.0	80	15	4	0.018"	C1725	 010764011955914817357503
357504	3.5	80	15	4	0.018"	C1725	 010764011955915517357504
357505	4.0	80	15	4	0.018"	C1725	 010764011955916217357505
357506	5.0	80	15	4	0.018"	C1725	 010764011955917917357506
366124	6.0	80	12	4	0.018"	C1725	 010764013041466217366124
366125	7.0	80	12	5	0.018"	C1725	 010764013041467917366125
357507	2.5	120	15	4	0.018"	C1725	 010764011955918617357507
357508	3.0	120	15	4	0.018"	C1725	 010764011955919317357508
357509	3.5	120	15	4	0.018"	C1725	 010764011955920917357509
357510	4.0	120	15	4	0.018"	C1725	 010764011955921617357510
357511	5.0	120	15	4	0.018"	C1725	 010764011955922317357511
366127	6.0	120	12	5	0.018"	C1725	 010764013041469317366127
366128	7.0	120	12	5	0.018"	C1725	 010764013041470917366128
366130	2.5	150	15	4	0.018"	C1725	 010764013041472317366130
366131	3.0	150	15	4	0.018"	C1725	 010764013041473017366131
366132	3.5	150	15	4	0.018"	C1725	 010764013041474717366132

Passeo-18 130 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
366133	4.0	150	15	4	0.018"	C1725	 010764013041475417366133
366134	5.0	150	12	4	0.018"	C1725	 010764013041476117366134
366135	6.0	150	12	5	0.018"	C1725	 010764013041477817366135
366136	7.0	150	12	5	0.018"	C1725	 010764013041478517366136
357512	2.5	170	15	4	0.018"	C1725	 010764011955923017357512
357513	3.0	170	15	4	0.018"	C1725	 010764011955924717357513
357514	3.5	170	15	4	0.018"	C1725	 010764011955925417357514
376292	4.0	170	13	4	0.018"	C1725	 010764013042999417376292
376293	5.0	170	13	4	0.018"	C1725	 010764013043000617376293
376294	6.0	170	12	5	0.018"	C1725	 010764013043001317376294
376295	7.0	170	12	5	0.018"	C1725	 010764013043002017376295
376297	2.5	200	14	4	0.018"	C1725	 010764013043004417376297
376298	3.0	200	14	4	0.018"	C1725	 010764013043005117376298
376299	3.5	200	14	4	0.018"	C1725	 010764013043006817376299
376300	4.0	200	13	4	0.018"	C1725	 010764013043007517376300
376301	5.0	200	13	4	0.018"	C1725	 010764013043008217376301
376302	6.0	200	12	5	0.018"	C1725	 010764013043009917376302
376303	7.0	200	12	5	0.018"	C1725	 010764013043010517376303

## Passeo-18 150 cm Catheter Length




















ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
366115	2.0	20	15	4	0.018"	C1725	 010764013041457017366115
366118	2.0	40	15	4	0.018"	C1725	 010764013041460017366118
366119	2.0	60	15	4	0.018"	C1725	 010764013041461717366119
366123	2.0	80	15	4	0.018"	C1725	 010764013041465517366123
366126	2.0	120	15	4	0.018"	C1725	 010764013041468617366126
366129	2.0	150	15	4	0.018"	C1725	 010764013041471617366129
366137	2.0	170	15	4	0.018"	C1725	 010764013041479217366137
376296	2.0	200	14	4	0.018"	C1725	 010764013043003717376296

## Passeo-35 Xeo Peripheral PTA Balloon Catheter










### Passeo-35 Xeo 90 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428777	3.0	20	21	5	0.035"	C1725	 010764013045148317428777
428778	4.0	20	20	5	0.035"	C1725	 010764013045149017428778
428779	5.0	20	19	5	0.035"	C1725	 010764013045150617428779
428780	6.0	20	16	5	0.035"	C1725	 010764013045151317428780
428781	7.0	20	15	5	0.035"	C1725	 010764013045152017428781
428782	8.0	20	14	6	0.035"	C1725	 010764013045153717428782









## Passeo-35 Xeo 90 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428783	9.0	20	12	6	0.035"	C1725	 010764013045154417428783
428784	10.0	20	11	6	0.035"	C1725	 010764013045155117428784
428786	3.0	40	21	5	0.035"	C1725	 010764013045157517428786
428787	4.0	40	20	5	0.035"	C1725	 010764013045158217428787
428788	5.0	40	19	5	0.035"	C1725	 010764013045159917428788
428789	6.0	40	16	5	0.035"	C1725	 010764013045160517428789
428790	7.0	40	15	5	0.035"	C1725	 010764013045161217428790
428791	8.0	40	14	6	0.035"	C1725	 010764013045162917428791
428792	9.0	40	12	6	0.035"	C1725	 010764013045163617428792
428793	10.0	40	11	6	0.035"	C1725	 010764013045164317428793
428794	12.0	40	10	7	0.035"	C1725	 010764013045165017428794
428795	3.0	60	21	5	0.035"	C1725	 010764013045166717428795
428796	4.0	60	20	5	0.035"	C1725	 010764013045167417428796
428797	5.0	60	19	5	0.035"	C1725	 010764013045168117428797
428798	6.0	60	16	5	0.035"	C1725	 010764013045169817428798
428799	7.0	60	15	5	0.035"	C1725	 010764013045170417428799
428800	8.0	60	14	6	0.035"	C1725	 010764013045171117428800
428801	9.0	60	12	6	0.035"	C1725	 010764013045172817428801
428802	10.0	60	11	6	0.035"	C1725	 010764013045173517428802




## Passeo-35 Xeo 90 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428803	12.0	60	10	7	0.035"	C1725	 010764013045174217428803
428804	3.0	80	21	5	0.035"	C1725	 010764013045175917428804
428805	4.0	80	20	5	0.035"	C1725	 010764013045176617428805
428806	5.0	80	19	5	0.035"	C1725	 010764013045177317428806
428807	6.0	80	16	5	0.035"	C1725	 010764013045178017428807
428808	7.0	80	15	5	0.035"	C1725	 010764013045179717428808
428809	8.0	80	14	6	0.035"	C1725	 010764013045180317428809
428810	9.0	80	12	6	0.035"	C1725	 010764013045181017428810
428811	10.0	80	11	6	0.035"	C1725	 010764013045182717428811
428812	12.0	80	10	7	0.035"	C1725	 010764013045183417428812
428814	4.0	100	20	5	0.035"	C1725	 010764013045185817428814
428815	5.0	100	19	5	0.035"	C1725	 010764013045186517428815
428816	6.0	100	16	5	0.035"	C1725	 010764013045187217428816
428817	7.0	100	15	5	0.035"	C1725	 010764013045188917428817
428818	8.0	100	14	6	0.035"	C1725	 010764013045189617428818
428819	9.0	100	12	6	0.035"	C1725	 010764013045190217428819
428820	10.0	100	11	6	0.035"	C1725	 010764013045191917428820
428821	12.0	100	10	7	0.035"	C1725	 010764013045192617428821
428823	4.0	120	20	5	0.035"	C1725	 010764013045194017428823




















### Passeo-35 Xeo 90 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428824	5.0	120	19	5	0.035"	C1725	 010764013045195717428824
428825	6.0	120	16	5	0.035"	C1725	 010764013045196417428825
428833	5.0	150	19	5	0.035"	C1725	 010764013045204617428833
428834	6.0	150	16	5	0.035"	C1725	 010764013045205317428834
428843	5.0	200	19	5	0.035"	C1725	 010764013045214517428843
428844	6.0	200	16	5	0.035"	C1725	 010764013045215217428844
428848	5.0	250	19	5	0.035"	C1725	 010764013045219017428848
428849	6.0	250	16	5	0.035"	C1725	 010764013045220617428849

















### Passeo-35 Xeo 130 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428851	3.0	20	21	5	0.035"	C1725	 010764013045222017428851
428852	4.0	20	20	5	0.035"	C1725	 010764013045223717428852
428853	5.0	20	19	5	0.035"	C1725	 010764013045224417428853
428854	6.0	20	16	5	0.035"	C1725	 010764013045225117428854
428855	7.0	20	15	5	0.035"	C1725	 0107640130452268 17428855
428856	8.0	20	14	6	0.035"	C1725	 010764013045227517428856
428857	9.0	20	12	6	0.035"	C1725	 010764013045228217428857
428858	10.0	20	11	6	0.035"	C1725	 010764013045229917428858

**Passeo-35 Xeo 130 cm Catheter Length (Continued)**

















ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428860	3.0	40	21	5	0.035"	C1725	 010764013045231217428860
428861	4.0	40	20	5	0.035"	C1725	 010764013045232917428861
428862	5.0	40	19	5	0.035"	C1725	 010764013045233617428862
428863	6.0	40	16	5	0.035"	C1725	 010764013045234317428863
428864	7.0	40	15	5	0.035"	C1725	 010764013045235017428864
428865	8.0	40	14	6	0.035"	C1725	 010764013045236717428865
428866	9.0	40	12	6	0.035"	C1725	 010764013045237417428866
428867	10.0	40	11	6	0.035"	C1725	 010764013045238117428867
428868	12.0	40	10	7	0.035"	C1725	 010764013045239817428868
428869	3.0	60	21	5	0.035"	C1725	 010764013045240417428869
428870	4.0	60	20	5	0.035"	C1725	 010764013045241117428870
428871	5.0	60	19	5	0.035"	C1725	 010764013045242817428871
428872	6.0	60	16	5	0.035"	C1725	 010764013045243517428872
428873	7.0	60	15	5	0.035"	C1725	 010764013045244217428873
428874	8.0	60	14	6	0.035"	C1725	 010764013045245917428874
428875	9.0	60	12	6	0.035"	C1725	 010764013045246617428875
428876	10.0	60	11	6	0.035"	C1725	 010764013045247317428876
428877	12.0	60	10	7	0.035"	C1725	 010764013045248017428877
428878	3.0	80	21	5	0.035"	C1725	 010764013045249717428878

## Passeo-35 Xeo 130 cm Catheter Length (Continued)









ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428879	4.0	80	20	5	0.035"	C1725	 010764013045250317428879
428880	5.0	80	19	5	0.035"	C1725	 010764013045251017428880
428881	6.0	80	16	5	0.035"	C1725	 010764013045252717428881
428882	7.0	80	15	5	0.035"	C1725	 010764013045253417428882
428883	8.0	80	14	6	0.035"	C1725	 010764013045254117428883
428884	9.0	80	12	6	0.035"	C1725	 010764013045255817428884
428885	10.0	80	11	6	0.035"	C1725	 010764013045256517428885
428886	12.0	80	10	7	0.035"	C1725	 010764013045257217428886
428887	3.0	100	21	5	0.035"	C1725	 010764013045258917428887
428888	4.0	100	20	5	0.035"	C1725	 010764013045259617428888
428889	5.0	100	19	5	0.035"	C1725	 010764013045260217428889
428890	6.0	100	16	5	0.035"	C1725	 010764013045261917428890
428891	7.0	100	15	5	0.035"	C1725	 010764013045262617428891
428892	8.0	100	14	6	0.035"	C1725	 010764013045263317428892
428893	9.0	100	12	6	0.035"	C1725	 010764013045264017428893
428894	10.0	100	11	6	0.035"	C1725	 010764013045265717428894
428895	12.0	100	10	7	0.035"	C1725	 010764013045266417428895
428896	3.0	120	21	5	0.035"	C1725	 010764013045267117428896
428897	4.0	120	20	5	0.035"	C1725	 010764013045268817428897



Passeo-35 Xeo 130 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428898	5.0	120	19	5	0.035"	C1725	 010764013045269517428898
428899	6.0	120	16	5	0.035"	C1725	 010764013045270117428899
428900	7.0	120	15	5	0.035"	C1725	 010764013045271817428900
428901	8.0	120	14	6	0.035"	C1725	 010764013045272517428901
428902	9.0	120	12	6	0.035"	C1725	 010764013045272517428901
428903	10.0	120	11	6	0.035"	C1725	 010764013045274917428903
428904	12.0	120	10	7	0.035"	C1725	 010764013045275617428904
428905	3.0	150	21	5	0.035"	C1725	 010764013045276317428905
428906	4.0	150	20	5	0.035"	C1725	 010764013045277017428906
428907	5.0	150	19	5	0.035"	C1725	 010764013045278717428907
428908	6.0	150	16	5	0.035"	C1725	 010764013045279417428908
428909	7.0	150	15	5	0.035"	C1725	 010764013045280017428909
428910	3.0	170	21	5	0.035"	C1725	 010764013045281717428910
428911	4.0	170	20	5	0.035"	C1725	 010764013045282417428911
428912	5.0	170	19	5	0.035"	C1725	 010764013045283117428912
428913	6.0	170	16	5	0.035"	C1725	 010764013045284817428913
428914	7.0	170	15	5	0.035"	C1725	 010764013045285517428914
428915	3.0	200	21	5	0.035"	C1725	 010764013045286217428915
428916	4.0	200	20	5	0.035"	C1725	 010764013045287917428916


















### Passeo-35 Xeo 130 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428917	5.0	200	19	5	0.035"	C1725	 010764013045288617428917
428918	6.0	200	16	5	0.035"	C1725	 010764013045289317428918
428919	7.0	200	15	5	0.035"	C1725	 010764013045290917428919
428920	3.0	250	21	5	0.035"	C1725	 010764013045291617428920
428921	4.0	250	20	5	0.035"	C1725	 010764013045292317428921
428922	5.0	250	19	5	0.035"	C1725	 010764013045293017428922
428923	6.0	250	16	5	0.035"	C1725	 010764013045294717428923
428924	7.0	250	15	5	0.035"	C1725	 010764013045295417428924












### Passeo-35 Xeo 170 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428927	5.0	20	19	5	0.035"	C1725	 010764013045298517428927
428928	6.0	20	16	5	0.035"	C1725	 010764013045299217428928
428929	7.0	20	15	5	0.035"	C1725	 010764013045300517428929
428930	8.0	20	14	6	0.035"	C1725	 010764013045301217428930
428935	4.0	40	20	5	0.035"	C1725	 010764013045306717428935
428936	5.0	40	19	5	0.035"	C1725	 010764013045307417428936
428937	6.0	40	16	5	0.035"	C1725	 010764013045308117428937
428938	7.0	40	15	5	0.035"	C1725	 010764013045309817428938

Passeo-35 Xeo 170 cm Catheter Length (Continued)




ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428939	8.0	40	14	6	0.035"	C1725	 010764013045310417428939
428940	9.0	40	12	6	0.035"	C1725	 010764013045311117428940
428941	10.0	40	11	6	0.035"	C1725	 010764013045312817428941
428944	4.0	60	20	5	0.035"	C1725	 010764013045315917428944
428945	5.0	60	19	5	0.035"	C1725	 010764013045316617428945
428946	6.0	60	16	5	0.035"	C1725	 010764013045317317428946
428947	7.0	60	15	5	0.035"	C1725	 010764013045318017428947
428948	8.0	60	14	6	0.035"	C1725	 010764013045319717428948
428949	9.0	60	12	6	0.035"	C1725	 010764013045320317428949
428950	10.0	60	11	6	0.035"	C1725	 010764013045321017428950
428953	4.0	80	20	5	0.035"	C1725	 010764013045324117428953
428954	5.0	80	19	5	0.035"	C1725	 010764013045325817428954
428955	6.0	80	16	5	0.035"	C1725	 010764013045326517428955
428956	7.0	80	15	5	0.035"	C1725	 010764013045327217428956
428957	8.0	80	14	6	0.035"	C1725	 010764013045328917428957
428958	9.0	80	12	6	0.035"	C1725	 010764013045329617428958
428959	10.0	80	11	6	0.035"	C1725	 010764013045330217428959
428971	4.0	120	20	5	0.035"	C1725	 010764013045342517428971

## Passeo-35 Xeo 170 cm Catheter Length (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
428972	5.0	120	19	5	0.035"	C1725	 010764013045343217428972
428973	6.0	120	16	5	0.035"	C1725	 010764013045344917428973
428980	4.0	150	20	5	0.035"	C1725	 010764013045351717428980
428981	5.0	150	19	5	0.035"	C1725	 010764013045352417428981
428982	6.0	150	16	5	0.035"	C1725	 010764013045353117428982
428990	4.0	200	20	5	0.035"	C1725	 010764013045361617428990
428991	5.0	200	19	5	0.035"	C1725	 010764013045362317428991
428992	6.0	200	16	5	0.035"	C1725	 010764013045363017428992
428995	4.0	250	20	5	0.035"	C1725	 010764013045366117428995
428996	5.0	250	19	5	0.035"	C1725	 010764013045367817428996
428997	6.0	250	16	5	0.035"	C1725	 010764013045368517428997

## Passeo-35 HP Peripheral High Pressure PTA Balloon Catheter




### Passeo-35 HP 40 cm Catheter Length

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
399063	6.0	40	25	6	0.035"	C1725	 010539152578127317399063
399067	7.0	40	23	6	0.035"	C1725	 010539152578131017399067
399071	8.0	40	22	6	0.035"	C1725	 010539152578135817399071

## Passeo-35 HP 75 cm Catheter Length












ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
399078	4.0	20	27	6	0.035"	C1725	 010539152578142617399078
399080	5.0	20	27	6	0.035"	C1725	 010539152578144017399080
399083	6.0	20	25	6	0.035"	C1725	 010539152578147117399083
399087	7.0	20	23	6	0.035"	C1725	 010539152578151817399087
399091	8.0	20	22	6	0.035"	C1725	 010539152578155617399091
399077	3.0	40	27	6	0.035"	C1725	 010539152578141917399077
399079	4.0	40	27	6	0.035"	C1725	 010539152578143317399079
399081	5.0	40	27	6	0.035"	C1725	 010539152578145717399081
399084	6.0	40	25	6	0.035"	C1725	 010539152578148817399084
399088	7.0	40	23	6	0.035"	C1725	 010539152578152517399088
399092	8.0	40	22	6	0.035"	C1725	 010539152578156317399092
399095	9.0	40	20	7	0.035"	C1725	 010539152578159417399095
399096	10.0	40	20	7	0.035"	C1725	 010539152578160017399096
399097	12.0	40	20	8	0.035"	C1725	 010539152578161717399097
399082	5.0	60	27	6	0.035"	C1725	 010539152578146417399082
399085	6.0	60	25	6	0.035"	C1725	 010539152578149517399085
399089	7.0	60	23	6	0.035"	C1725	 010539152578153217399089
399093	8.0	60	22	6	0.035"	C1725	 010539152578157017399093

## Passeo-35 HP 75 cm Catheter Length (Continued)










ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
399094	8.0	80	22	6	0.035"	C1725	 010539152578158717399094
399086	6.0	100	25	6	0.035"	C1725	 010539152578150117399086
399090	7.0	100	23	6	0.035"	C1725	 010539152578154917399090


## Fortress Peripheral Reinforced Introducer Sheath System

### Fortress Introducer Sheath

ORDER NUMBER	FRENCH SIZE (F)	SHEATH LENGTH (CM)	SHEATH SHAPE	DILATOR(S) INCLUDED	UNITS PER BOX	C-CODE	BARCODE (GTIN)
358813	4.0	45	Curved	Curved/ Straight	5	C1894	 010425124450362817358813
358814	4.0	45	Straight	Curved/ Straight	5	C1894	 010425124450360417358814
444486	4.0	65	Straight	Straight	5	C1894	 010425124450363517444486
444485	4.0	90	Straight	Straight	5	C1894	 010425124450364217444485
386590	5.0	45	Curved	Straight	5	C1894	 010425124450367317386590
386591	5.0	45	Straight	Straight	5	C1894	 010425124450365917386591
444484	5.0	65	Straight	Straight	5	C1894	 010425124450057317444484
444483	5.0	90	Straight	Straight	5	C1894	 010425124450369717444483
386593	6.0	45	Curved	Straight	5	C1894	 010425124450372717386593
386594	6.0	45	Straight	Straight	5	C1894	 010425124450370317386594
444482	6.0	65	Straight	Straight	5	C1894	 010425124450373417444482

## Fortress Introducer Sheath (Continued)

ORDER NUMBER	BALLOON DIAMETER (MM)	BALLOON LENGTH (MM)	RATED BURST PRESSURE (ATM)	SHEATH COMPATIBILITY (F)	RECOMMENDED GUIDEWIRE (IN)	C-CODE	BARCODE (GTIN)
444481	6.0	90	Straight	Straight	5	C1894	 010425124450374117444481
452386	7.0	45	Curved	Straight	5	C1894	 010425124450377217452386
452382	7.0	45	Straight	Straight	5	C1894	 010425124450375817452382
452383	7.0	65	Straight	Straight	5	C1894	 010425124450378917452383
452388	7.0	90	Straight	Straight	5	C1894	 010425124450379617452388
452387	8.0	45	Curved	Straight	5	C1894	 010425124450382617452387
452385	8.0	45	Straight	Straight	5	C1894	 010425124450092417452385
452384	8.0	65	Straight	Straight	5	C1894	 010425124450095517452384
452389	8.0	90	Straight	Straight	5	C1894	 010425124450096217452389



# Innovation When It Matters Most

Orsiro Mission, PK Papyrus, Pantera, Oscar, Pulsar, Astron, Fortress and Passeo are trademarks or registered trademarks of the BIOTRONIK Group of Companies.

For Professional Use Only. Product details presented here are intended for healthcare professionals in the United States.

Please refer to each product's Instructions for Use on the BIOTRONIK manuals website ([manuals.biotronik.com](http://manuals.biotronik.com)) for a complete listing of the indications, contraindications, warnings and precautions.

BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97035-5369  
(800) 547-0394 (24-hour)  
(800) 291-0470 (fax)  
[www.biotronik.com](http://www.biotronik.com)  
©2024 BIOTRONIK, Inc.  
All rights reserved.  
11/13/24 VIP\_002\_24e

Rx  
ONLY



ReCross, NHancer, TrapIT and MicroRx are trademarks or registered trademarks of IMDS.  
Manufactured by IMDS and distributed by BIOTRONIK in selected countries.



**BIOTRONIK**  
excellence for life